

Upload an Omnipod[®] PDM to Insulet Provided Glooko with Glooko Kiosk

What You Need:



— or —



Omnipod DASH[™] System PDM

Omnipod[®] System PDM

Micro to Micro Cable or
Micro to Mini Cable

Glooko Kiosk

Follow These Steps:

1 Upload Omnipod[®] PDM Data

Omnipod DASH[™] System:



- Turn DASH[™] on, and unlock it
- Plug DASH[™] into the Kiosk with the Micro-to-Micro cable
- Tap **Export** on DASH[™]
- Press **Start Here** on the Kiosk
- Select **Omnipod DASH System**
- Follow the prompts, and press **Sync** on the Kiosk

NOTE: The patient will receive an email to complete activation of their account. That way, they will automatically be connected to your account so you can see their data when they upload from home.

Omnipod[®] System:



- Plug the Omnipod[®] System into the Kiosk with the Micro-to-Mini USB cable
- Press **Start Here** on the Kiosk
- Select **Omnipod System**
- Follow the prompts, and press **Sync** on the Kiosk

2 View Omnipod[®] PDM Data

- Log into your account on my.glooko.com to view the uploaded data

Still experiencing issues after attempting these instructions?

Please contact Glooko Support.

Call: 800.206.6601 ▪ Text: 650.720.5310 ▪ Email: support@glooko.com

For more information on the Omnipod[®] System, contact the Insulet Customer Care Team:
800.591.3455 | myomnipod.com

Insulet Corporation
100 Nagog Park
Acton, MA 01720
978.600.7850

All screen images are examples, and for instructional use only.