


Upload an Omnipod[®] PDM to Insulet Provided Glooko with the Glooko Transmitter

What You Need:



Follow These Steps:

1 Upload to Glooko Transmitter

- Ensure the Glooko Transmitter displays **Ready to Transfer**.

Omnipod DASH[™] System:




- Turn DASH[™] on, and unlock it
- Plug DASH[™] into Transmitter with the USB-A to Micro cable
- Tap **Export** on DASH[™]
- Unplug DASH[™], and plug it back in
- Once the upload is complete, the Glooko Transmitter will display **Transfer Complete**

NOTE: If user does not know the passcode, use the last 4 digits of the serial number found on the back of DASH[™].

Omnipod[®] System:



- Plug the Omnipod[®] System into the Transmitter with the USB-A to Mini cable
- The Omnipod[®] System will beep 4 times, then display **USB Device Ready**
- Once the upload is complete, the Glooko Transmitter will display **Transfer Complete**




- Log into your account on **my.glooko.com** to view the uploaded data

2 Assign Device and View Data on Computer

- Go to **my.glooko.com** and log into your Glooko account.
- In the Assign Devices tab, check the box in the right-hand column next to the device that was uploaded.

Note: devices will be displayed in the Assign Devices tab for 24 hours, so be sure to assign data to correct patient within this timeframe



(Step 2 continued on other side)

2 Assign Device and View Data on Computer (cont'd)

For Existing Patients:

1. VIEW MATCHES:

If the device has been uploaded before, you will see the patient name that matches the device serial number. Select the patient name and click **Assign**.

- Choose to either create a PDF report, or view patient's data online

2. SEARCH FOR A PATIENT:

If no matches were found automatically, use the search box to type the patient's name, then click **Assign**.

- Choose to either create a PDF report, or view patient's data online

For New Patients:

1. CREATE A NEW ACCOUNT:

Select **Create New Account** and fill out the patient's information.

- Choose to either create a PDF report, or view patient's data online

NOTE: The patient will receive an email to complete activation of their account. That way, they will automatically be connected to your account so you can see their data when they upload from home.

Still experiencing issues after attempting these instructions?

Please contact Glooko Support.

Call: 800.206.6601 ▪ Text: 650.720.5310 ▪ Email: support@glooko.com

For more information on the Omnipod® System, contact the Insulet Customer Care Team:
800.591.3455 | myomnipod.com

Insulet Corporation
100 Nagog Park
Acton, MA 01720
978.600.7850

All screen images are examples, and for instructional use only.