


THE LOGBOOK REPORT

Jen Levi

Date of Birth: 5/4/2001 Diabetes: Type 1

Logbook

May 1st - May 14th, 2017 (14 days)




The **Logbook Report** shows consolidated data including BG readings, insulin delivery, carb intake, and pod changes for a customized time frame.

Jen Levi

Date of Birth: 5/4/2001 Diabetes: Type 1

Logbook

1 May 1st - May 14th, 2017 (14 days)



Glucose

- Above Range (>180 mg/dL)
- In Target Range (70 -180 mg/dL)
- Below Range (<70 mg/dL)
- 123 Manual BG

Bolus

- Before Meal
- After Meal
- Multiple Readings (most extreme shown)

Bolus

- Override
- Multiple Boluses
- Extended Bolus
- Manual Bolus

Basal

- Basal Rate
- > Temp. Basal Rate
- ↔ Suspend
- Set/Site Change

3

Carb

FEATURES

- 1 **Customized time frame**
Example: May 1 - May 14 (14 days).
- 2 **Glucose averages**
Example: there were 42 BG readings recorded from May 1 - May 14, with an average reading of 124 mg/dL.
- 3 **Average total daily dose (TDD)**
Example: On average, 40 units of insulin were used per day from May 1 - May 14.
- 4 **Color-coded blood glucose readings**
Example: a reading that is above target range (182) was recorded around 11a.m. on Tuesday, May 14.
- 5 **Customized time frame, with most recent day at the top**
Example: Pod change on Mon 5/13.

THE SUMMARY REPORT


Larry Jones

Date of Birth: 1/1/1977 Diabetes: Type 1


Summary

May 6th - June 4th, 2017 (30 days)


GLUCOSE (CGM)



INSULIN



Bi-Hourly



DIET

Carbs/Day 123 g
 Entries/Day 2.6

FITNESS

Steps/Day 3,729


COMMENTS

The comments will be displayed in this area

Time of Day



Day of Week



The **Summary Report** displays glucose readings that are above, below and in range, as well as the insulin delivery for the customized time frame.

Larry Jones

Date of Birth: 1/1/1977 Diabetes: Type 1


Summary

1 May 6th - June 4th, 2017 (30 days)

GLUCOSE (CGM)





INSULIN



DIET
Carbs/Day: 123 g
Entries/Day: 2.6

FITNESS
Steps/Day: 3,729

COMMENTS
The comments will be displayed in this area



FEATURES

- 1 Customized time frame**
Example: May 6 - June 4 (30 days).
- 2 Glucose averages and statistics**
Example: 6% of CGM readings between May 6 - June 4 are below range, with the lowest reading recorded as 32 mg/dL.
- 3 Insulin usage**
Example: On average, the daily insulin delivery breakdown between May 6 - June 4 is:
 - Basal: 58% (21.8 units)
 - Bolus: 42% (15.8 units)
 - Total Daily Dose (TDD): 37.6 units
- 4 Glucose averages displayed in three different time segments (bi-hourly, time of day, and day of week)**
Example: In the bi-hourly time segment view, the average CGM reading recorded between 12 a.m. and 2 a.m. from May 6 - June 4 is 163.


THE OVERVIEW REPORT

Jenn Lee

Date of Birth: 5/4/2001 Diabetes: Type 1

Overview

May 28th - Jun 10th, 2017 (14 days)



GLUCOSE (CGM)					GLUCOSE (BG)		INSULIN			OTHER	
126 mg/dL	72 mg/dL	3%	87%	10%	125 mg/dL	36 mg/dL	21.5 units	55% (11.9 U)	45% (9.6 U)	61,865	118 g
Average	Standard Deviation	Low < 70 mg/dL	In Target Range 70-180 mg/dL	High > 180 mg/dL	Average	Standard Deviation	Total Insulin/Day	Total Bolus/Day	Basal/Day	Exercise (steps)	Carbs/Day

Glucose

- Before Meal Target Range (70-130 mg/dL)
- After Meal Target Range (70-180 mg/dL)
- Above Range
- In Target Range
- Below Range


Types of Readings

- ● BG Readings
- ● Pump BG
- ◆ ◆ Calibrations
- CGM Readings
- ◆ Pump BG Above 400

Bolus

- Bolus
- Basal

The **Overview Report** shows consolidated daily data including BG and CGM readings, carb intake, insulin delivery and exercise (if applicable).



FEATURES

- 1 Customized time frame**
Example: May 28 - June 10 (14 days).
- 2 BG readings and CGM calibrations, overlaid with CGM trend line**
- 3 Glucose averages and statistics**
Example: 3% of CGM readings from May 28 - June 10 are below target range.
- 4 Insulin usage**
Example: On average, 21.5 units of insulin were delivered per day from May 28 - June 10.
- 5 Additional information**
Note: steps available if user connects a compatible activity tracker
Example: On average, 118 grams of carbs were entered per day from May 28 - June 10.


THE OVERLAY REPORT— Spaghetti and AGP *(for use with CGM)*

Larry Jones

Date of Birth: 1/1/1977 Diabetes: Type 1

Overlay - Spaghetti & AGP

May 1st - 14th, 2017 (14 days)



GLUCOSE (CGM)

169 mg/dL Average	70.6% % Time CGM Active	72 mg/dL Standard Deviation	46.3% Coefficient of Variation	9% Very Low < 54 mg/dL	10.1% Low < 70 mg/dL	64.6% In Target Range 70-180 mg/dL	25.1% High > 180 mg/dL	6.8% Very High > 250 mg/dL
----------------------	-------------------------------	--------------------------------	-----------------------------------	------------------------------	----------------------------	--	------------------------------	----------------------------------

Glucose

- Target Range (70-180 mg/dL)
- Monday
- Tuesday
- Wednesday
- Thursday
- Friday
- Saturday
- Sunday

Percentiles

- 25-75%
- 10-90%
- Median
- Lowest - Highest

The Overlay – Spaghetti & AGP Report shows:

Spaghetti Graph (*top*): View CGM trend lines, color coded by each day of a customized time frame


Ambulatory Glucose Profile (AGP) Report (*bottom*): View CGM readings from multiple days overlaid in a single, 24-hour view.

Larry Jones

Date of Birth: 1/1/1977 Diabetes: Type 1

Overlay - Spaghetti & AGP

1 May 1st - 14th, 2017 (14 days)



GLUCOSE (CGM)		
169 mg/dL Average	70.6% % Time CGM Active	72 mg/dL Standard Deviation
	46.3% Coefficient of Variation	

9% Very Low < 54 mg/dL	10.1% Low < 70 mg/dL	64.6% In Target Range 70-180 mg/dL	25.1% High > 180 mg/dL	6.8% Very High > 250 mg/dL
------------------------------	----------------------------	--	------------------------------	----------------------------------

Glucose
 — Target Range (70-180 mg/dL)
 --- Thursday
 --- Monday
 --- Tuesday
 --- Wednesday
 --- Friday
 --- Saturday
 --- Sunday

Percentiles
 ■ 25-75%
 ■ 10-90%
 — Median
 --- Lowest - Highest

FEATURES

1 Customized time frame

Example: Report includes data from May 1 - May 14 (14 days).

2 Glucose averages and statistics

Example: The standard deviation based on the CGM data from May 1 - May 14 is 72 mg/dL.

3 Percent of time spent low, high and in range


Example: 9% of CGM readings from May 1 - May 14 were very low (lower than 54 mg/dL).

4 AGP Report Key

- **Dark blue line:** median curve, showing median glucose values at that time point
- **Darker shade:** 25th-75th percentile (IQR) showing half of all glucose readings
- **Lighter shade:** 10th-90th percentile which captures glucose excursions
- **Dotted lines:** lowest and highest readings at that time of day

Example: Based on the 14 days from May 1 - May 14, most readings around 3 p.m. were above target range.

THE OVERLAY REPORT— Scatter & Box Plot (for SMBG)




The Overlay Scatter & Box Plot Report shows manual BG readings plotted by time of day, and color coded based on target range.

Larry Jones

Date of Birth: 1/1/1977 Diabetes: Type 1

Overlay - Scatter & Box Plot

1 May 1st - 14th, 2017 (14 days)



GLUCOSE (BG)							0.8	
169 mg/dL	72 mg/dL	46.3%	9%	10.1%	64.6%	25.1%	6.8%	BG Tests/Day
Average	Standard Deviation	Coefficient of Variation	Very Low < 54 mg/dL	Low < 70 mg/dL	In Target Range 70-180 mg/dL	High > 180 mg/dL	Very High > 250 mg/dL	

2

Glucose
 - - - Before Meal Target Range (70-130 mg/dL) ● Above Range
 — After Meal Target Range (70-180 mg/dL) — Median
 ● In Target Range

FEATURES

1 Customized time frame

Example: Report includes data from May 1 - May 14 (14 days).

2 Glucose averages and statistics

Example: The average BG readings from the selected 14 days is 169.


THE CALENDAR REPORT

Larry Jones


Date of Birth: 1/1/1977 Diabetes: Type 1

Calendar

May 1st - 14th, 2017 (14 days)



The Calendar Report shows consolidated daily data including BG and CGM readings and calibrations, carb intake, insulin delivery and pod changes.



FEATURES

1 Customized time frame

Example: Report includes data from May 1 - May 14 (14 days).

2 Day of the week

Example: data in this box is from May 1.

3 Bolus and basal unit scale

4 BG readings color-coded based on target range

Example: a BG reading above target range (303) was recorded on Wednesday, May 3.

5 CGM calibrations

Example: a CGM calibration was recorded the morning of Sunday, May 7.

6 Basal profile

Example: This line shows the basal profile from May 1 - May 7.

7 Pod change

Example: There was a pod change on Thursday, May 4 around 3 p.m.

How to Print and Store a PDF Report


CREATE A REPORT:

From the Population Tracker:

1. Right click on a patient's name
2. Select **Print Report**
3. Choose desired reports and time frame
4. Click **Create**

From Patients Profile:

1. Click **Create Report** in the top right-hand corner
2. Choose desired reports and time frame
3. Click **Create**



STORE A REPORT:

Once the report is created, download it to your desktop. Then, either:

1. Upload the report to your EHR system.
- or
2. Print and scan the report to your EHR system.