



Omnipod[®] 5 & Glooko[®] Key information

AN OMNIPOD[®] 5 LEAFLET FOR
PEOPLE LIVING WITH DIABETES



LINK your Omnipod 5 Automated Insulin Delivery System data with Glooko

If you have missed linking your Omnipod 5 and Glooko, please follow these four simple steps:

1. Log into your online Omnipod account at www.omnipod.com, using your Omnipod ID
2. Go to the **My Account** tab and then select Account linking from the drop down menu. Then click “**Account Linking**”
3. Select “**sign up**” or “**log in**” with Glooko
4. Proceed to Glooko and complete the process to link your account

When complete, your status will update to “**Linked**” and data will automatically upload into your personal Glooko account*. You can always unlink your account between Omnipod and Glooko to stop data sharing.



NOTE: All screen images are examples, for instructional purposes only.

*Must be connected to WiFi or cellular data, and Omnipod 5 users must have WiFi or cellular data access to receive updates in Glooko.

Choose to **SHARE** your
Omnipod 5 Automated
Insulin Delivery
System data with your
healthcare provider

If you have missed linking
your Glooko account to your
Healthcare Professional, please
follow these five simple steps:



1. Log into your Glooko account at my.glooko.com
2. Select “**Settings**” from the drop-down menu (top right on the home screen)
3. Scroll down to Account and select “**+ Add New Code**”
4. Select “**Continue**”
5. Enter ProConnect Code (listed below) then select “**Submit**”


*Entering a ProConnect code in
Glooko is the only way for your
healthcare provider to see your
Omnipod 5 data on their computer.*

**ProConnect Code:
Provided by your
healthcare practitioner**

NOTE: All screen images are examples, for instructional purposes only.

Using the GLOOKO MOBILE APP

*You have the choice to use the
Glooko mobile app or the online
Glooko platform my.glooko.com*



NOTE: All screen images are examples, for instructional purposes only.

For iOS



Download the app and login

- From the App Store search for and install the “Glooko mobile app”
- Log in or create an account by filling out the requested information


View your data

Once your data is uploaded to Glooko, you can view your trends and graphs on your iPhone or iPad

Choose to share data with your healthcare provider

- Click the **More** tab,
- Click **Share Reports** to email or fax the reports

For Android



Download the app and login

- From the Google Play™ Store search for and install the “Glooko mobile app”
- Log in or create an account by filling out the requested information

Sync and view your data

- Tap Sync from the Home Screen, and select Omnipod 5 Automated Insulin Delivery System
- Follow the prompts to complete the sync

Choose to share data with your healthcare provider

Click **Share Reports** to email or fax the reports.

NOTE: All screen images are examples, for instructional purposes only.

Create, Store or Print a **PDF REPORT**

You can create and view reports on the **Glooko mobile app** or at **my.glooko.com**.

1. Log in to your account

2. To create a report:

- Click on your **Profile**
- Click **Create PDF Report** (top right)
- Choose desired reports and time frame and click **Create PDF**

3. To save/print a report:

Once the report is created, either save it to your desktop or a desired location, or print directly.

Create PDF Close

Jane Doe Preferred PDF Setting
DOB: 17/08/09 No Profile Available
Diabetes: Type 1 Manage Favourites

Time: 2 weeks 10/11/2021 - 23/11/2021

Summary Logbook Overview Daily Overview Week View

Overlay 2 pages Calendar 2 pages Insights 1 page Devices 14 pages

Estimated report length: 20

Type comment here! Comments will appear in the Summary report section.

Print mode: Black and White Colour

Save selection as Favourite Profile Create PDF

NOTE: All screen images are examples, for instructional purposes only.

TIP! Save Customised Report Sets:


- Click **Settings**
- Scroll down and click **New Favourite**
- Select desired reports, and name the report set (*i.e. Dr. Smith reports*)



Great reports to get started with are:


Summary reports, Daily, Weekly

Snapshot of **key statistics** and **aggregated information** to identify positives and challenges to support your care



The Device Settings Report

When you upload your data to Insulet Provided Glooko, your pump settings will be saved. This report is a great resource for you and your care team to reference.



NOTE: All screen images are examples, for instructional purposes only.



If you have any questions or concerns regarding Omnipod 5 please contact the Omnipod team on:

0800 011 6132*

or +44 20 3887 1709 if calling from abroad.

If you are experiencing any issues with Glooko, please contact Glooko Support:

020 7795 8191*

www.glooko.com
uk@glooko.com



*Your call be monitored and recorded for quality monitoring purposes. Calls to 0800 numbers are free from local landlines, but other networks may charge for these calls.

©2024 Insulet Corporation. Omnipod and the Omnipod logo are trademarks or registered trademarks of Insulet Corporation in the United States of America and other various jurisdictions. Glooko is a trademark of Glooko, Inc. and used with permission. All rights reserved. [Insulet Netherlands BV Stadsplateau 7, 3521 AZ Utrecht, Netherlands]. INS-OHS-11-2023-00141 V2