



INSTRUCTIONS FOR USE
GLOOKO® FOR PERSONAL USE

Contents

Intended use.....	3
Intended users	3
Contraindications.....	3
Clinical benefits	3
Warnings.....	3
1. How to Start Using Glooko®	4
1.1 Performance characteristics	4
1.2. Compatibility	4
Operating Systems	4
Device compatibility.....	5
IT and security	5
MeterSync Blue*	5
Glooko Uploader	5
1.3. Upload and View Your Diabetes Data	6
Step 1: Log into your Glooko Account	6
Step 2: Upload your Diabetes Device Data	7
Step 3: View your Diabetes Data.....	10
2. Glooko Mobile App Overview	11
2.1. Manage your Account.....	12
Profile.....	12
Add a ProConnect Code.....	12
Devices.....	13
Connect to Apple Health	13
Settings	14
Prime Detection*	14
2.2. Manually add data	15
2.3. Set Reminders	16
2.4. View Charts and Graphs.....	16
Graphs.....	17
Trends*	17
2.5. Synchronise Fitness & Health Apps	18
3. Glooko Web App Overview.....	19
3.1. Manage Settings.....	19
My Profile	19
Account.....	20

Change Email, Password or Language	20
Add a ProConnect Code.....	20
Change Your Primary Provider	21
Communications	21
Data sharing & privacy	21
Data Settings	22
Apps & Devices.....	23
Connect Apps & Devices.....	23
Upload Omnipod® PDM	24
Connect to DASH™ Cloud*	24
Get Glooko Uploader	24
Connect Eversense Account.....	25
Connect iGlucose Devices	25
Connect Abbott FreeStyle LibreView account*	26
Upload LibreView data*	26
Account Management	27
Delete Account.....	27
4. Reports Overview	28
4.1. Available Reports.....	28
4.2. Create and Share Reports	29
Mobile App Reports	29
Web App Reports	29
5. Support	30
Appendix: Additional Features	31
A.1. Glooko Clinical Research*	31
A.2. Basal-IQ.....	31
A.3. Control-IQ	32
A.4. Omnipod 5 System*	32
A.5. CamAPS FX*	33

Intended use

The Glooko Mobile Application is a diabetes management platform that can be used at home by patients. It allows viewing of patient data from compatible medical devices and other health monitoring devices. It also allows Health Care Professionals to support the treatment of diabetes through analysis and monitoring of a patient's diabetes therapy and Care Program. Glooko Mobile is intended for use by individuals with diabetes.

Intended users

The Glooko Mobile Application is intended for use by individuals with diabetes. The intended patient population is individuals with diabetes.

Contraindications

None known.

Clinical benefits

Patients who use the Glooko System may experience the following clinical benefits:

- Improvement in glycaemic control
- Improved A1C
- Sustained glycaemic improvements

Warnings

Glooko does not measure, interpret or make decisions on the data it conveys, nor is it intended to provide automated treatment decisions or be used as a substitute for professional judgement. All medical diagnosis and treatment are to be performed under the supervision and oversight of an appropriate healthcare provider.

If you have concerns about your current diabetes diagnosis and treatment, please contact your healthcare provider.

1. How to Start Using Glooko®

Glooko allows people with diabetes to synchronise data from popular blood glucose (BG) meters, insulin pumps and Smart Pens, continuous glucose monitors (CGMs) and health and fitness devices – and securely share that information with their healthcare teams. With multiple viewing options, Glooko provides insights into glucose data fluctuations based on integrated or manually added lifestyle events, including exercise, food, carbs and medications. Glooko can be accessed via the Glooko mobile app or Glooko web app.

NOTE: Some screens may vary based on additional features. For more information, consult the [Appendix: Additional Features](#).

1.1 Performance characteristics

The Glooko Web Application maintains a minimum service uptime of 99%.

1.2. Compatibility

Operating Systems

For each Operating System and web browser below, Glooko supports the latest version as well as the version preceding that (as of this document's release date).

Mobile Operating Systems:

- iOS
- Android

Computer Operating Systems:

- macOS
- Windows

Web browsers:

- Google Chrome
- Microsoft Edge
- Firefox
- Safari

For the best performance, we recommend using Google Chrome.

Device compatibility

To determine your diabetes device compatibility, consult the [Glooko Compatibility](#) page.

The complete Glooko System is comprised of the following Glooko products:

- Glooko Web Application (REF-0001)
- Glooko Mobile Application (REF-0002)
- Glooko Research Application (REF-0009)

**May not be available in your country.*

IT and security

Glooko requires a username and password for accessing a Glooko Web and Mobile account. Do not share your username or password with anyone.

Glooko Web and Mobile are connected to the Internet (Wi-Fi or cellular data network).

MeterSync Blue*

The MeterSync Blue is a Bluetooth device used with BG meters that are not Bluetooth-enabled to upload diabetes device data to the Glooko mobile app. For detailed instructions on how to pair the MeterSync Blue to your mobile device, consult the [MeterSync Blue Quick Start Guide](#).

** Currently only available in the U.S.*

Glooko Uploader

If your diabetes device is not compatible with your smartphone, and the PC/MAC option is available according to the Compatibility page, you can install the Glooko Uploader to upload your data using your computer. For detailed instructions on how to use this tool to upload data, consult the [Glooko® Uploader - Instructions for Use](#).

1.3. Upload and View Your Diabetes Data

Once you have determined your diabetes device compatibility and which method you will use to upload your data, follow the steps below to begin uploading your diabetes data to Glooko:

Step 1: Log into your Glooko Account



Log into your Glooko account via the Glooko mobile app or Glooko web app.

Mobile App: Install the app from the [App Store](#) or [Google Play](#). If you already have an account, tap **Log In**, enter your login information and tap **Log In** again. If you do not have an account, tap **Register** and follow the on-screen prompts to create a new account. For regulatory compliance, minors will need parental or guardian consent to activate an account.

NOTE: If you were given an activation code by your provider or through a sponsored programme, you will have the option to enter that code during the registration process.

Web App: Go to my.glooko.com in your web browser. If you already have an account, enter your login information and click **Sign In**. If you do not have an account, click the link to sign up for a new account (or go to join.glooko.com) and follow the on-screen prompts to create a new account. For regulatory compliance, minors will need parental or guardian consent to activate an account.

NOTE: If you received an activation request via email notifying you that your healthcare provider created an account for you, click the **Activate Account** button within the email and follow the on-screen prompts to activate your account. If you were given an activation code by your provider, go to join.glooko.com/activate and follow the on-screen prompts to activate your account.



Step 2: Upload your Diabetes Device Data

After you have logged into the Glooko mobile or web app, you can upload your data to your Glooko account.


NOTE: Having the wrong time or date on your diabetes device can cause issues when trying to sync with Glooko. When first viewing your data within the Glooko Mobile or Web App, confirm that the data timestamps within the app correspond to the timestamps on your diabetes device. In the **Devices** tab in the Glooko Web App, you can confirm that there is no time offset alert (clock icon) for physically uploaded devices.

The screenshot shows the Glooko web application interface. At the top, the 'glooko' logo is on the left, and a user profile icon 'LT' is on the right. Below the header is a navigation bar with tabs: Summary, Graphs, Logbook, History, Insights, and **Devices** (which is highlighted with an orange box). The main content area displays the profile for 'Lucy Test' with details like DOB (01/01/1970) and Diabetes (Type 1). Below the profile are links for 'Create PDF Report' and 'Upload Omnipod® PDM'. A blue 'Add device' button is on the right. Under 'My cloud connections', a 'Dexcom' connection is shown as 'Connected' with a last sync time of 31 Oct 2022 - 10:52. Under 'My devices', an 'Ascensia CONTOUR NEXT ONE' device is listed with ID 7802HP211215 and a last sync time of 23 Feb 2023 - 16:15. To the right of this device entry, a larger box shows the device name 'Ascensia CONTOUR NEXT ONE' next to a clock icon and a time offset '+00:59', which is highlighted with an orange box. Below this, it states 'Device settings are not available'.

Mobile app:

The Glooko Mobile App has been redesigned and users will now be encouraged to complete a set of onboarding goals in order to quickly familiarise themselves with the app and its main features.

NOTE: The onboarding goals may vary depending on account type and country/language settings.



To add and sync additional devices after the onboarding goals have been completed, tap the **Sync device** button on the home screen. Please reference the [Glooko Mobile App Overview](#).

Enable automatic reminders for glucose testing and syncing:*

Users who have a BG meter connected, but are not using a CGM, will be invited to the **Smart reminders** feature via a pop-up dialogue within the Glooko Mobile app. Please follow the on-screen instructions to enable reminders.

** Currently only available in the US.*

Order hardware:*


If you need to order syncing hardware, tap the **Account** menu in the upper-right corner and then select **Order Sync Hardware**.

**Only applicable in the U.S. and Canada*

Web App:

Some diabetes devices can be synced using your computer together with the [Glooko Uploader](#) software to upload data to your Glooko account, while others require a cloud connection or can be

synced using a mobile phone and the Glooko Mobile App. To find out the different uploading methods available for each diabetes device, please click the **Add device** button on the Devices page (or use the Add device link on the Summary page).



1. Click the **Add device** button on the Devices page.
2. Select device category.
3. Select your device and upload data according to the available methods listed for the device.
4. The devices that you have previously synced with your account will be visible under My devices or My cloud connections depending on device.


Please note: The compatibility of the diabetes devices varies depending on country/market.

Step 3: View your Diabetes Data

Once your data has been uploaded, it is available for you to view in the Glooko mobile app or Glooko web app.

Toggle between the different screens for multiple views of your data – and generate [reports](#) to print, share or save that data.

To learn more, see [View Charts and Graphs](#).



NOTE: In the Glooko Web App, insulin pump data and Smart Pen data will be displayed as **Insulin – Device** and manually entered insulin data will be displayed as **Insulin - Manual**. If a patient has a pump with advanced technology to automatically stop the pump and/or automatically deliver insulin, the summary page will display an additional information card with the title **System details** below Insulin. Please reference [Appendix: Additional Features](#) for more information.

2. Glooko Mobile App Overview

There are a number of Glooko mobile app features that add value to your diabetes management programme and enable Glooko to more seamlessly integrate with your lifestyle. By taking advantage of these features, you can easily track, upload and share your diabetes data wherever you are.

Add care team

Tap here to add the ProConnect Code of your care team to share your data with them.

Daily log

Tap here to see daily historic data, such as: Blood Glucose Readings, Insulin Doses, Food Entries, Exercise data and more.


Get Started*

Complete the onboarding goals by tapping each card and following the on-screen instructions. Once all the goals have been completed this section will not appear again.

**The onboarding goals may vary depending on account type and country/language settings.*

Add

Tap here to make manual entries, such as: A1C, Blood Pressure, Exercise, Food, Insulin, Medication, Note and Weight.



Account

Tap here to find the following menu options: Profile, Devices, Settings, Share report, Reminders, Communications, Data sharing & privacy, Order sync hardware, About and Help.

Sync device

Tap here to add and sync your diabetes devices.

Charts and graphs

Tap here to view charts and graphs.

BG/CGM, insulin, exercise, food, weight and blood pressure

If data exists from within the last two weeks, a short summary will be presented here for these data types.

View previous weeks

Tap here to see a weekly overview, including: Time in Range/week, average Blood Glucose/week, average Insulin/day.

Top trend*

A recent significant blood glucose trend is displayed here. Tap **View all** to see more trends.

**Currently only available in the US.*

2.1. Manage your Account

By tapping the **Account** icon at the top-right corner of the home screen you can view and update your personal details, add a ProConnect Code, manage your device settings, share reports, create reminders, delete or log out of your account and more.

NOTE: Many of these and additional settings can be configured in the Glooko web app > [Settings](#).

Profile

To access your personal details, select **Profile** from the **Account** menu.

To update your information, tap an item to change the selected value. Tap **Done**, if applicable, to register your changes.

NOTE: Name, Email and Gender can only be updated in the Glooko web app > [Settings](#).

Add a ProConnect Code


If your healthcare provider(s) also use Glooko, you can add their ProConnect Code(s) to your account to automatically share your diabetes data with them via the secure cloud.

NOTE: You can add up to five (5) ProConnect Codes to your account.

To add a new ProConnect Code:

1. In **Account** > **Profile**, tap **ProConnect**.
2. Tap **Add Care Team**.
3. Enter the ProConnect Code and then tap **Continue**
4. Tap **Connect**.
5. Tap **OK** to save.

TIP: You can also add a ProConnect Code in the Glooko web app in Settings > [Account](#).




The screenshot shows the 'Profile' screen in the Glooko app. At the top, there is a back arrow on the left and a close 'X' icon on the right. Below the title bar, the profile information is listed in a table-like format with labels on the left and values on the right. The fields are: Name (Lucy Test), Email (blue-thompson-360@example.com), Type of Diabetes (Type 1), Gender (Female), Date of Birth (1 January 1970), Height (167.6 cm), Weight (72.6 kg), and Country (United States of America). Below these fields, there are two rows with labels 'ProConnect' and 'Delete Account', each followed by a right-pointing chevron. At the bottom of the screen, there is a 'LOG OUT' button.

Profile	
Name	Lucy Test
Email	blue-thompson-360@example.com
Type of Diabetes	Type 1
Gender	Female
Date of Birth	1 January 1970
Height	167.6 cm
Weight	72.6 kg
Country	United States of America
ProConnect	>
Delete Account	>
LOG OUT	

Devices

In **Account > Devices**, you can set up a [MeterSync Blue](#) device, view Meter Time Settings, connect an Apple Health account and toggle MeterSync cable syncing on or off.

To update your Devices, tap an item and change the settings as needed. Tap **Done**, if applicable, to register your changes.



Connect to Apple Health

You can connect your Apple Health account to your Glooko account to synchronise BG and CGM readings, insulin data and nutrition information in near real-time.

NOTE: Before connecting your Apple Health account, the Apple Health app must be installed on your mobile device. You will be required to make updates within the Apple Health app to allow Glooko to access your data.

To connect your Apple Health account:

1. In **Account > Devices**, tap **Apple Health Settings**.
2. Tap on the Shared Data Types (**Carbs, Calories, Fat, Protein, Insulin** or **Blood Glucose**) you would like to enable.

NOTE: Tapping **Blood Glucose** will enable sharing of BG and CGM data.


3. When the Apple Health app opens:
 - Tap **Turn All Categories On** or tap to toggle on synchronisation for only **Blood Glucose, Carbohydrates, Dietary Energy, Insulin Delivery, Protein** or **Total Fat**.
 - Tap **Allow** to allow Glooko to access your data.
4. You will be routed back to the Glooko app to view your synchronised data.

Settings

In **Account > Settings**, you can view and update your Glooko account settings, including your Glucose Data Source, Prime Detection, Unit of Measurement*, Set Daily Time Ranges and Set Target BG Range. You can also enable data sharing with Dexcom.

To update your settings, tap an item to change the selected value. Tap the **tick (✓)** or select **Done**, if applicable, to register your changes.


** For users in Canada, this setting is locked to mmol/L and cannot be toggled to mg/dL.*



Prime Detection*

For other Smart Pens

This is a feature available for patients using selected Smart Pens. Prime detection is automatically enabled for patients where this feature is allowed by local regulations in the patients' country of residence. When prime detection is enabled, Glooko excludes probable prime doses from the insulin statistics and graphs. All doses are however stored in Glooko and are never deleted. In the **Daily log**, patients have the option to manually change the status of a dose.



NOTE: The priming detection algorithm is not always accurate. The software detects prime doses that are two units or less within six minutes before another insulin injection. If too much time has passed between your priming dose and injected dose, your prime dose may show up as an injected dose. You will have to manually designate your dose as prime. To do this, tap on the dose in the **Daily log** and select **Primed**.

For Sanofi SoloSmart/Novo Nordisk Mallya

The device detects prime doses and reports them to Glooko. For detailed information regarding which doses the SoloSmart/Mallya reports as prime doses, please refer to the SoloSmart/Mallya instruction manual.

** May not be available in your country. Currently not available in the U.S.*

2.2. Manually add data


Different types of data, including A1c, Blood Pressure, Weight, Medication, Insulin, Food and Exercise, as well as custom Notes, can be manually added to your Glooko mobile app so that they can be viewed alongside blood glucose data to reveal trends and patterns.

To manually add data:

1. Tap **Add (+)** on the home screen.
2. Select the type of data: **A1c, Blood Pressure, Weight, Medication, Insulin, Food, Exercise or Notes**.
3. Tap the **date/time** field and adjust the date and time of the data entry, if necessary.
4. Search for and/or select (from your Recent, Frequent or Custom lists) the item(s) you would like to include in the entry and add all relevant details. Tap **Done** and/or the **plus (+)** symbol to add the item(s).
 - To add multiple entries (of different types) at the same time, toggle between the event types at the top of the Add Event menu and add items as needed (not applicable for A1c, Blood Pressure or Weight).
5. Tap the **tick (✓)** at the top-right of the screen to save your entries.

TIP: When adding a Food event, you can search for the item, select from Recent or Frequent items, add a Custom item or add a quick Carbs entry. You can also use your voice or the bar-code scanner by tapping the **microphone** icon (🎤) or **bar-code** icon (📊) in the search bar. As foods are added, total carbs will display at the bottom of the screen (expand for calories, fat and protein).

NOTE: The ability to search for food items is currently available in English, German, French, Spanish and Italian. And the microphone and bar code search options are currently available in the U.S. only. In other languages you have the option to add a Custom item or a Carbs entry.



2.3. Set Reminders

You can set reminders to Check Your Blood Sugar, Take Insulin, Take Medication or set a Custom reminder from within the Glooko mobile app. You will be alerted via push notification on your mobile device at the specified reminder time.

NOTE: In order to receive push notifications from Glooko, you must also toggle on the Push Notifications option in Profile > [Settings](#).

To create a reminder:

1. Select **Reminders** from the **Account** menu.
2. Tap the **plus symbol (+)** at the top-right of the screen to create a new reminder.
3. Configure the reminder settings:
 - **Time of Day:** Tap and scroll to select the time of day the reminder should occur.
 - **Day(s) of the Week:** Tap to place a tick (✓) beside the day(s) of week on which the reminder should occur.
 - **Select Reminder Type(s):** Tap to place a tick (✓) beside the reminder type(s) you would like to add and enter all relevant information.
4. Tap **Save** at the top-right of the screen to enable the reminder(s).

The screenshot shows the 'Set Reminders' interface. At the top, there is a back arrow, the title 'Set Reminders', and a 'SAVE' button. Below the title, there are two sections: 'Schedule Time' and 'Set Reminders'. The 'Schedule Time' section has a dropdown menu showing '11:35 pm' and another dropdown showing 'Everyday'. The 'Set Reminders' section contains a list of reminder types, each with a checkbox and an icon. The types are: 'Test your Blood Sugar' (green icon, checked), 'Take Insulin' (blue icon, checked), 'Insulin' (blue icon, unchecked, with a dropdown for 'units'), 'Take Medication' (purple icon, unchecked), and 'Custom' (blue icon, checked). At the bottom, there is a field labeled 'Enter custom reminder' with a plus sign to its right.

2.4. View Charts and Graphs

Glooko makes visualising your diabetes health easier by providing at-a-glance insights into your blood glucose data in the form of Graphs and Trends – and you can generate [reports](#) to save, print or share that data.

You can access this view by tapping **Charts and Graphs** on the home screen.


Graphs

The Graphs tab displays visual representations of your BG/CGM, Exercise, Insulin and Carb data. Graphs can be viewed over a 1-week, 2-week, 1-month or 3-month time period.

To access Graphs:

1. Tap **Charts and Graphs** on the home screen.
2. Tap the **Graphs** tab at the top of the screen.
3. The following graphs display:
 - BG Average or CGM Average
 - Blood pressure
 - Exercise
 - Insulin
 - Carbs
 - BG by Time of Day or CGM Daily Overview
 - Weight
4. Tap the **down arrow (▼)** to expand any of the graphs, and slide to the left or right to view greater detail by day.

TIP: If you have BG and CGM data available, you can select which Glucose Data Source will pull in Profile > [Settings](#).




Trends*

The Trends tab displays Patterns, Comparisons and Best Days based on your BG/CGM data. Trends can be viewed over a 1-week, 2-week, 1-month or 3-month time period.

NOTE: A minimum of 15 days of data, with at least two readings per day, is required to provide insights.

To access Trends:

1. Tap **Charts and Graphs** on the home screen.
2. Tap the **Trends** tab at the top of the screen.
3. The following Trends display:
 - **Patterns:** These identify times of day when your BG/CGM readings are notably high, low or in your target range. Patterns are available only in the 1-month date range and will populate only if they can be identified.



- **Comparisons:** These help you monitor changes in your diabetes data over time by indicating how often your BG/CGM readings were high, low or in range for the time period selected as compared with the previous time period.
 - **Best Days:** These help you identify which days your diabetes health was best managed, based on your BG/CGM data.
4. Tap the **down arrow (▼)** to expand a Comparison to view greater detail.

TIP: If you have BG and CGM data available, you can select which Glucose Data Source will pull in Profile > [Settings](#).

** Currently only available in the U.S.*

2.5. Synchronise Fitness & Health Apps

With the Glooko mobile app, you have the ability to integrate exercise routine data from leading third-party fitness and health apps. This allows you to visualise your activity data and how it correlates with your glucose and insulin data.

To find out if your fitness or health app(s) are compatible with Glooko, consult the [Glooko Compatibility](#) page.

To synchronise your fitness or health app(s) with Glooko:

1. Tap the **Devices** button on the home screen.
2. Tap **Add New Device**.
3. Tap to place a **tick (✓)** beside **Fitness and Health Devices**, then tap **Continue**.
4. Tap to place a **tick (✓)** beside the device(s) you would like to synchronise.
5. Tap **Add**.
6. Next, tap **Connect** to complete the setup process for each new device.
 - After adding your device(s), you can complete the setup process at a later time by tapping **Set Up** beside the device(s) on the Sync menu.
7. Follow the on-screen prompts to connect the device(s) to Glooko.

3. Glooko Web App Overview

The Glooko web app enables access to your glucose, insulin, carb and fitness data through consolidated analytics and report views. With graphs and statistics by time of day and date range, you can better understand how often your glucose readings fall below, above and within normal range as a result of your lifestyle and care decisions. Reports are based on data from your blood glucose meter, insulin pump, continuous glucose monitor (CGM), fitness and health apps and events manually added in the Glooko mobile app.

3.1. Manage Settings

To access your account Settings, select **Settings** from the drop-down menu at the top-right of your Glooko web app. On this screen, you have the option to update your Glooko account information (including your Email, Password and ProConnect Codes), configure your Data Settings, connect to third-party apps and more.

The screenshot shows the Glooko web app interface. At the top, there's a blue header with the 'glooko' logo on the left and a user profile icon labeled 'LT' on the right. Below the header is a navigation bar with links: Summary, Graphs, Logbook, Insights, and Devices. A dropdown menu is open from the user icon, showing options: Lucy Test, Settings (highlighted with an orange box), Help, and Log Out. The main content area is titled 'My Profile' and contains a form with the following fields: First Name (Lucy), Surname (Test), Date of Birth (January 1, 1970), Type of Diabetes (Type 2), Gender (Female), Height (168 cm, with radio buttons for feet and cm), and Weight (73 kgs, with radio buttons for lbs and kgs). A blue 'Save' button is located at the bottom right of the form, also highlighted with an orange box.

My Profile

In Settings > My Profile, you can view and update basic demographic information, including your Name, Type of Diabetes, Date of Birth, Height and Weight. To make changes, click into the field(s) and update as needed. Click **Save** to register your changes.

Account

In Settings > Account, you can update your Email, Password or Language, manage your ProConnect Code(s) and change your Primary Provider.

The screenshot shows the 'Account' settings page. At the top, there are three sections: 'Email Address' with the value 'blue-thompson-360@example.com' and a 'Change Email' link; 'Password' with masked characters '*****' and a 'Change Password' link; and 'Language' with a 'Change Language' link. Below these is the 'ProConnect Code' section, which includes a question mark icon and an 'Edit ProConnect' link. The 'Primary Provider' section is marked with a star and a question mark icon, showing 'Diabetes Clinic' as the selected provider with a code '1234567' and a lock icon. There is a '+ Add New Code' button below the provider list. The 'Data Authorization' section at the bottom contains two rows of radio buttons: 'Share Omnipod® System data with Insulet' with 'Agree' and 'Disagree' options, and 'Share Glooko data with Dexcom' with 'Agree' (selected) and 'Disagree' options. Each row has a question mark icon for help.

Change Email, Password or Language

To update your Email, Password or Language, click **Change Email**, **Change Password** or **Change Language**, enter or select the desired value and click **Save**.

Add a ProConnect Code

If your healthcare provider(s) also use Glooko, you can add their ProConnect Code(s) to your account to automatically share your diabetes data with them via the secure cloud.

NOTE: You can add up to five (5) ProConnect Codes to your account.

To add a new ProConnect Code, follow these steps:

1. In Settings > Account, click **+ Add New Code**.
2. A pop-up window will display. Click **Continue** to proceed.
3. Enter the ProConnect Code.
4. Click **Submit**.
5. Click **OK** to return to the Settings menu.

TIP: You can also add a ProConnect Code to your account in the Glooko mobile app in Profile > [Personal](#).

Change Your Primary Provider

Your Primary Provider is indicated with a star (★). Your Primary Provider has the ability to adjust your Pump BG Entry Settings, which influence your graphs and statistics. If you are sharing with multiple providers, only one can be set as your Primary Provider.

To change which provider is designated as your Primary Provider, follow these steps:

1. In Settings > Account, click **Edit ProConnect**.
2. Select your choice for **Primary Provider** by toggling on the button to the left of the provider name.
3. Click **+ Confirm Change** to register the new selection.

Communications

In **Settings > Communications**, you have the option to manage how you receive messages and updates from the app. Under **Mobile phone experience**, you can enable or disable SMS notifications to receive reminders and messages tailored to your app usage and connected devices. In the **Marketing** section, you can choose whether you want to stay updated on the latest developments, new features and company news. Please note that essential account-related emails, such as password reset notifications, will still be sent even if you opt out of marketing updates.

Data sharing & privacy

In **Settings > Data sharing & privacy**, you can manage how your data is shared and review privacy policies. Under **Device data sharing**, you can choose to share your data with Glooko's partners to help improve products and services, troubleshoot devices or contribute to health research projects. If you agree, the selected partners will access and use your data in accordance with their privacy policies. Additionally, you can opt into the **Glooko Patient Research Network** to learn about potential research opportunities that match your profile. Finally, you can review Glooko's **Terms of Use** and **Privacy Policy** for detailed information about data usage.

Data Settings

In Settings > Data Settings, you have the option to adjust your Unit of Measurement, Pump BG Entry Settings, Target BG Ranges and Daily Time Ranges.

The following Settings display:

- **Unit of Measurement:** This can be toggled to **mg/dL** or **mmol/L**. This is automatically set based upon whether you are using a mg/dL or mmol/L-based BG meter. For users in Canada, this setting is locked to mmol/L and cannot be toggled to mg/dL.

- **Pump BG Entry Settings:** This can be toggled to **Yes** or **No** to include or exclude BG readings manually entered in insulin pumps in your graphs and statistics. By default, these readings are included.

NOTE: If you are [ProConnected](#), your Primary Provider can also adjust your Pump BG Entry Settings.

- **Target BG Ranges:** This allows you to change your Target BG Ranges, including the Lower Limit, Before Meal Upper Limit and After Meal Upper Limit. To edit these ranges, click into the field of the value that needs to be modified and enter the desired value.
- **Daily Time Ranges:** This allows you to set the times of day that indicate the start of your Morning, Afternoon, Evening and Night routines. To update these ranges, click the **down arrow (▼)** beside any of the times of day and select a new start time.


TIP: Many of these settings can also be updated in the Glooko mobile app in Profile > [Settings](#).

The screenshot shows the 'Data Settings' window with the following sections:

- Unit of Measurement:** Radio buttons for 'mg/dL' (selected) and 'mmol/L'.
- Pump BG Entry Settings:** Radio buttons for 'Include in statistics' (selected 'Yes') and 'No'. A 'Restore Default' link is on the right.
- Target BG Ranges:** Input fields for 'Lower Limit' (70 mg/dL), 'Before Meal Upper Limit' (130 mg/dL), and 'After Meal Upper Limit' (180 mg/dL). A 'Restore Default' link is on the right.
- Daily Time Ranges*:** A sequence of colored arrows representing the day: 5:00 AM (Morning), 10:00 AM (Afternoon), 3:00 PM (Evening), 9:00 PM (Night), and Morning. Each has a dropdown arrow. A 'Restore Default' link is on the right.
- Footer:** A small note: '*Events that occur before Morning will appear in the previous Night. To confine events strictly to calendar days, set Morning to start at midnight.' and a blue 'Save' button.

Apps & Devices

In **Settings > Apps & Devices**, you have the option to connect your third-party fitness and health apps, install the Glooko Uploader or connect a device account to your Glooko account. Please note that the ability to connect different apps and device accounts will vary depending on region/market. You can also manage your apps and devices in the **Devices** tab.



Connect Apps & Devices

You have the ability to synchronise your third-party fitness and health app(s) with Glooko to pull in exercise routine data.

To connect your third-party app(s) to Glooko:

1. In Settings > Apps & Devices, click **Connect Apps & Devices**.
2. Choose from the list of third-party apps and click **Connect**. You will be redirected to the login page for the app(s).
3. Follow the on-screen prompts to log in and connect the app(s) to your Glooko account.

Upload Omnipod® PDM

To upload data from your Omnipod system:

1. In Settings > Apps & Devices, click **Upload Omnipod® PDM**.
2. When prompted to confirm if you would like to proceed, click **Continue**.
3. Select **Omnipod® PDM System** or **Omnipod DASH™ System**, then click **Next**.

NOTE: If you are attempting to synchronise an Omnipod DASH™ System on a Mac computer, you will be prompted to install the [Glooko Uploader](#) to synchronise your data to Glooko.

4. Follow the on-screen prompts to connect your device and upload your data.

TIP: If you have synchronised an Omnipod device previously, you have the option to upload data by selecting the **Upload Omnipod® PDM** option at the top-right of most screens in your account. You can also synchronise your Omnipod device in the Glooko mobile app via Home > [Sync](#).

Connect to DASH™ Cloud*

You can connect your Omnipod DASH™ Cloud account to your Glooko account to synchronise your data once every 24 hours.

To connect your Omnipod DASH™ Cloud account:

1. In Settings > Apps & Devices, click **Connect to DASH™ Cloud**.
2. Click **Connect**.
 - If you have an existing PodderCentral™ account, enter your login information and click **Log In**.
 - If you do not have an existing PodderCentral™ account, scroll down to **Create a PodderCentral™ Account** and create a new account.
3. Follow the on-screen prompts to register and/or connect your account.

NOTE: After connecting your account, the device status will display as Pending, and it may take up to an hour for data to appear. Once the status changes to Connected, your data will begin synchronising.

TIP: You can also connect your DASH™ Cloud account in the Glooko mobile app via Home > [Sync](#).

**May not be available in your country.*

Get Glooko Uploader

You can install the Glooko Uploader to synchronise data from compatible blood glucose meters, insulin pumps and CGMs to your Glooko account from your computer. To learn more, click [here](#).

Connect Eversense Account

You can connect your Eversense account to your Glooko account to synchronise data from Eversense once every four hours.

To connect your Eversense account:

1. In Settings > Apps & Devices, click **Connect Eversense Account**.
2. Enter your Eversense Email and Password, then click **Sign In**.
 - If you do not have an account, click **Sign Up** to create a new account.
3. Follow the on-screen prompts to create and/or connect your account.

TIP: You can also connect your Eversense account in the Glooko mobile app via Home > [Sync](#).

**May not be available in your country.*

Connect iGlucose Devices

You can connect your iGlucose meter(s) to your Glooko account to synchronise data from iGlucose in real time.

NOTE: Before connecting your iGlucose meter to Glooko, you must have at least one reading on your meter.

To connect your iGlucose meter:

1. In Settings > Apps & Devices, click **Connect iGlucose Devices**.
2. Enter your iGlucose meter's Serial Number and Last Reading (value), then click **Next**.
3. Follow the on-screen prompts to connect your meter.

TIP: You can also connect your iGlucose meter in the Glooko mobile app via Home > [Sync](#).

Connect Abbott FreeStyle LibreView account*

You can connect your Abbott FreeStyle Libreview account to your Glooko account to continuously sync data from your Abbott FreeStyle Libre sensor.

To connect your Abbott FreeStyle Libreview account:

1. In Settings > Apps & Devices, click **Connect Abbott FreeStyle Libreview account**.
2. Enter your email address and your FreeStyle Libre app Password, then click **Sign In**.
 - If you do not have an account, go to <https://www.libreview.com> to create a new account.
3. Follow the on-screen prompts to create and/or connect your account.

**Currently available only in Sweden and Norway.*

Upload LibreView data*

You can add data from your Abbott FreeStyle Libre device to Glooko by downloading your Abbott data in LibreView as a CSV file and then uploading the CSV file to your Glooko account.

To upload a CSV file with LibreView data to your Glooko account:

1. In **Settings > Apps & Devices**, click **Upload LibreView data**.
2. Follow the on-screen prompts, then click **Upload**.


** May not be available in your country.*

Account Management

Delete Account

In Settings > Account Management, you have the option to delete your Glooko account at any time. This will also delete all of your Glooko account data. To delete your account, simply click the **Delete Account** button.

TIP: You can also delete your account in the Glooko mobile app via Profile > [Personal](#).



4. Reports Overview

With Glooko, it's easy to view and share your diabetes data. If there is data available in your Glooko account, you can save, print and share reports via the Glooko mobile app or Glooko web app.

For detailed information about the available reports and how to interpret the data, view the [Glooko Report Reference Guide](#).

4.1. Available Reports

Available reports include:

- Summary (only available in the Glooko web app)
- Logbook
- Overview
- Daily Overview
- Overlay
- Week View (only available in the Glooko web app)
- Calendar
- Insights (only available in the Glooko web app)
- Devices

4.2. Create and Share Reports


When creating reports, you can include all available data or select which data should display.

Mobile App Reports

In the Glooko mobile app, you have the option to Email, Fax or Print a PDF report or Email your report in CSV format.

To create reports in the Glooko mobile app:

1. Select **Share Report** from the side navigation bar (≡) on Android devices or the More menu on iOS devices.
2. Tap to select your desired sharing method.
 - If sharing a PDF, select the timeframe, your report criteria and preferred print option.
 - Tap the **Reports** (📄) icon at the top-right of the screen
3. Follow the on-screen prompts to generate your report.




Web App Reports

In the Glooko web app, you have the option to save or print a report in PDF format.

To view and share reports in the Glooko web app:

1. Select **Create PDF Report** at the top-right of most screens.
2. Select a time range and which report data should be included in the report.
3. Once the report criteria are defined, click **Create PDF**.
4. View, download or print your report.



5. Support

If you have questions, we're always happy to help. Our Support Team is available Monday through Friday, 8 a.m. until 8 p.m. ET. You can reach out to us in any of the following ways:

- Web Support: support.glooko.com
- Email Support: support@glooko.com
- SMS Support: +1-650-720-5310

Any serious incident that has occurred in relation to the device should be reported to the Glooko Support Team and the competent authority of the country in which you are established.

In the event where the Glooko System is malfunctioning, please discontinue use of the product and contact Glooko Support. Examples of malfunctions include incorrect device time data and incorrect units of measurement.



Glooko Inc.
579 University Avenue
Palo Alto, CA 94301
USA



Glooko AB
Nellickevägen 20
SE-412 63 Göteborg
Sweden



REF-0001



REF-0002



REF-0009

Appendix: Additional Features


A.1. Glooko Clinical Research*

Glooko Clinical Research allows healthcare providers to monitor users participating in clinical research studies. As a research participant, you have access to a separate Glooko Clinical Research mobile app that must be connected to a unique Glooko Clinical Research ProConnect Code (provided by your site coordinator) in order to share your diabetes data.

As a research participant you can also upload and share your data, using a separate Research Uploader.

NOTE: The Glooko Clinical Research mobile app can be connected to only one ProConnect Code at a time.

**May not be available in your country.*




A.2. Basal-IQ

Patients using a Tandem t:slim X2 pump together with a CGM have the option to activate a technology called Basal-IQ. It is an advanced technology that predicts and helps prevent low blood sugar. Basal-IQ allows a patient to let the pump automatically suspend and resume insulin delivery based on the CGM readings.

A patient with Basal-IQ installed on the pump can choose between 2 modes:

- Basal-IQ
- Manual Mode

In Glooko, if Basal-IQ data exists for a patient, this is featured on the Summary page in form of an information card called: **System Details**. Basal-IQ is also presented in the day view of the **Graphs** section.




A.3. Control-IQ

Patients using a Tandem t:slim X2 pump together with a CGM have the option to activate a technology called Control-IQ. It is an advanced hybrid closed-loop technology that predicts and helps prevent both highs and lows. Control-IQ allows a patient to let the pump automatically adjust insulin levels based on the CGM readings.

A patient with Control-IQ installed on the pump can choose between 4 modes:

- Control-IQ
- Sleep
- Exercise
- Manual




In the Glooko web app, if Control-IQ data exists for a patient, this is featured on the Summary page in form of an information card called: **System Details**. Control-IQ is also presented in the day view of the **Graphs** section.

A.4. Omnipod 5 System*

Patients using an Omnipod 5 pump together with a Dexcom CGM have the option to activate the closed-loop technology developed by Insulet. This is an advanced hybrid closed-loop technology that predicts and helps prevent both highs and lows. It allows a patient to let the pump automatically adjust insulin levels based on the CGM readings.

A patient with an Omnipod 5 and a Dexcom CGM can choose between 4 modes:

- Automated mode
- Automated: Limited
- Automated: Activity
- Manual mode



In Glooko, if closed-loop data from the Omnipod 5 System exists for a patient, this is featured on the Summary page in form of an information card called: **System Details**. Closed-loop data is also presented in the day view of the **Graphs** section.

NOTE: The Omnipod 5 System is a cloud-to-cloud integration. The connection with Glooko is established at www.omnipod.com, where the users need to authorize the connection between their Omnipod 5 System and Glooko. Once the connection is made, data will stream regularly into Glooko with an hour delay.


**May not be available in your country.*

A.5. CamAPS FX*

Patients using a pump compatible with both Glooko and the CamAPS FX app have the option to use CamAPS FX's advanced adaptive hybrid closed-loop technology, which automatically adjusts insulin delivery to the insulin pump based on the sensor glucose readings.

The following 5 modes are visualised in Glooko:

- Auto mode 'On'
- Boost
- Ease-off
- Attempting
- Auto mode 'Off'



In Glooko, if closed-loop data from the CamAPS FX app exists for a patient, this is featured on the Summary page in form of an information card called: **System details**. Closed-loop data is also presented in the day view of the **Graphs** section.

**May not be available in your country.*