

GET STARTED with Omnipod[®] 5



STEPS TO COMPLETE BEFORE ATTENDING OMNIPOD 5 TRAINING

- Do you have a vial of rapid-acting U-100 insulin ready for training?
- Are you either ready to wear/wearing a Dexcom sensor (with your Dexcom receiver switched off) or prepared to attend with a new unopened FreeStyle Libre 2 Plus sensor ready for activation during training?
- Do you have your Omnipod[®] ID and password ready for training? Follow the steps in this guide for online Omnipod[®] Account creation
- Bring a fully charged Controller with you when attending training

Important Considerations

COMPLETE OMNIPOD 5 ONBOARDING


This guide will provide support for creating an Omnipod 5 and completing the necessary consent screens to get started on Omnipod 5.

SENSOR REQUIREMENTS

Using Dexcom

If you plan to use Omnipod 5 with a Dexcom sensor* you will need to attend the training already wearing a Dexcom sensor which must be used with a corresponding compatible Dexcom mobile app.

Note: Dexcom receivers must be switched off prior to training.



Using FreeStyle Libre

When using the Omnipod 5 system with the FreeStyle Libre 2 Plus sensor*, all sensor-specific actions and alerts are controlled through your Omnipod 5 Controller. You cannot use another device to search, read values, or receive alarms from the sensor.

Note: If attending an Omnipod 5 training make sure to attend the training session with a new unopened FreeStyle Libre 2 Plus sensor.



CONNECT TO OMNIPOD DISCOVER™

Omnipod Discover is the diabetes data management system that connects to Omnipod 5. By connecting your Omnipod 5 account with Omnipod Discover you will provide permission and enable data sharing† allowing you to:

- Review your personalised diabetes insights including insulin delivery insights, trends and patterns
- Share data with your Healthcare professional and/or diabetes clinic
- Receive weekly detailed reports from Discover providing overall insulin usage and readings from your glucose monitoring sensor

Recommendation: Set up the Controller's Wi-Fi to automatically connect to your home or work network. This will help ensure it connects regularly—ideally as often as possible, but at least once per week.†

Screen image for illustrative purposes only. Pod and sensors shown without necessary adhesive.

*Compatible sensors are sold and prescribed separately. Sensor availability may vary by country or region.

†Connections between Omnipod 5 Controller and Insulet Cloud are via Wi-Fi.



How to Get Started on Omnipod 5

Before attending Omnipod 5 training you will need to create your online Omnipod Account.

WHY SHOULD YOU REGISTER?

An online Omnipod Account is required to use Omnipod 5. Once you have created your online Omnipod Account and provided the necessary consents you will be able to review orders, track shipments, and request support if needed.

During Onboarding you will have the opportunity to connect your Omnipod Account with Omnipod Discover. Providing access to additional diabetes data and insights and enabling you to share your Omnipod 5 system data with your Healthcare team.

HOW TO REGISTER

Visit **www.omnipod.com/setup** or scan this QR Code from a smartphone.

In your Omnipod 5 Starter Kit you will receive a Welcome Card containing the same QR code and website address for account set up. To get started and create an online Omnipod Account you will need an e-mail address.



WHEN SHOULD YOU REGISTER?

Online Omnipod Account registration must be completed before attending Omnipod training when you will activate your Omnipod 5 Controller for the first time.

1 - Starting Omnipod Onboarding

Go to **www.omnipod.com/setup** or scan the QR code.
Then, select your country.



To login for the first time, click on "Create Omnipod ID" and follow the steps.


Confirm the Region of Residence

18+ YEARS

UNDER 18 YEARS

Note: Make sure the e-mail address provided is correct. If there are errors or typos in the email address used you will not be able to complete Omnipod 5 Onboarding.

Complete the fields. If a parent or legal guardian is completing on behalf of a dependent, you will need to create the Omnipod ID for yourself using your own details.



Click on the link
in the e-mail.




Review and confirm information.
Make sure the correct country code is
showing in the phone number field.

2 - Creating an Omnipod ID and Password

Create Omnipod ID and password



The screenshot shows a web form titled "Create your personal Omnipod® ID". It includes fields for "First name", "Last name", "Email", and "Password". To the right of the form, there are two sections: "Omnipod® ID requirements" and "Password requirements", each with a list of rules and green checkmarks indicating compliance. At the bottom of the form are "Back" and "Next" buttons.

Omnipod ID requirements:

- Must be 6 characters long
- Contains no special characters
- Contains no blank spaces

Password requirements:

- Must be 8 characters long
- Contains no blank spaces
- Must include an uppercase, lowercase, number and one of the following special characters !#\$%+<->@_
- Cannot include part of your first name, last name or Omnipod ID/Username


Select your method of two-factor authentication:

- by E-mail (Required)
- by SMS (Optional)



The screenshot shows a web form titled "Sign in security". It has two radio button options: "Email (Required)" and "Text Message (SMS) (Optional)". At the bottom are "Back" and "Next" buttons.

The Omnipod ID and password have been created! Now it's time to finalise registration.





Enter the Omnipod ID and password.



If this is the first time logging in, or if logging in from a different device, follow the steps to verify the account.

Please note: If this is entered incorrectly after 5 attempts, you will be locked out of your account for 30 minutes. Try again after 30 minutes has passed.

Once the verification code is received as per your chosen method. Enter this code into the field as prompted to proceed.



3 - Reading and Validating Consents

Data Privacy

At Insulet, the safety and security of our users and products is paramount in everything we do. We are dedicated to making the lives of people with diabetes easier and simplify diabetes management. Insulet respects the privacy of every one of our patients and is committed to the protection of their personal information. We have dedicated teams that are focused on keeping patient information safe from unauthorized access.

During the following steps, a number of data privacy consents will need to be validated: conditions of use of the product, the management of personal data, the use of the device and the improvement of products and services.

Omnipod® 5 Terms & Conditions

In order to start using your Omnipod® 5 Automated Insulin Delivery System you must acknowledge that you have read and agree to Insulet Corporation's Terms and Conditions which include the following:

[Go to User License Agreement, Omnipod 5 System Warning, Terms of Use](#)

Agree



Omnipod® 5 Consents
Data Privacy Introduction

Thank you for choosing the Omnipod® 5 Automated Insulin Delivery System. We respect the importance of your privacy. As part of your Omnipod® 5 onboarding, we therefore provide you on the following page with information on the processing of your personal data and ask for you to respond to various consent requests for such processing. We refer to our Privacy Policy for further information.

You can change your choices at any time in your online Omnipod® account by contacting the Omnipod® Customer Care Team at 1-800-522-2147 or +1-212-522-2147 toll-free.

*Your call may be recorded for quality monitoring and training purposes. Calls to 800-numbers are free from local landlines but other networks or network carriers may charge for these calls.

Continue



Omnipod® 5 Consents
Product Use

I agree that Insulet may collect and process my personal data, including my health data, as explained to me in the Insulet Privacy Policy for the purpose of the provision and use of the Omnipod® 5 Automated Insulin Delivery System, including verification and authentication of my Omnipod® 5 and activation of the Omnipod® 5 System, the provision of product support and maintenance of the Omnipod® 5 System.

Please contact the Omnipod® Customer Care Team if you have any questions at 1-800-522-2147 or +1-212-522-2147 toll-free.

*Your call may be recorded for quality monitoring and training purposes. Calls to 800-numbers are free from local landlines but other networks or network carriers may charge for these calls.

Agree

If the Product Research, Development and Improvement consent is accepted a list of optional questions will be shown.

Omnipod® 5 Consents
Product Research, Development and Improvement

I agree that Insulet may collect and process my personal data, including my health data, to research and improve the Omnipod® 5, Automated Insulin Delivery System, and related products and services, and to help Insulet's development of new products and services, as explained to me in the Insulet Privacy Policy. I understand that this may also include combining my data with other data Insulet holds on me, such as my demographic data, and data Insulet receives from my continuous glucose monitor, provider and/or diabetes data management provider.

Note: This consent is optional and you can use the Omnipod® 5 System without providing such consent. Not providing this consent will not impact your ability to use the Omnipod® 5 System. If you consent to such processing, you will be asked in the next step to optionally provide additional personal data.

Agree and Continue
Skip and Continue

By answering a few quick questions about your diabetes history and therapy, you'll help us better understand how Omnipod 5 works for different people—whether newly diagnosed or managing diabetes for years. Your input supports important research, improvements to support to your needs, and contributes to the development of future products and services.

Omnipod® 5 Consents
Customer Health Data

Gender
☐ Female
☐ Male

Are you currently using Omnipod®?
☐ No
☐ Yes

What type of diabetes therapy were you using before you chose to go on Omnipod® 5?

What type of glucose sensor were you using before using Omnipod® 5?

What type of diabetes do you have?

How long have you had diabetes?

What was your last HbA1c result?


Agree and Continue
Skip and Continue



4 - Omnipod Discover and Account Linking

Omnipod Discover is a diabetes data management system with a simple user interface designed specifically for people with diabetes. Omnipod 5 users can track glucose and insulin delivery patterns to highlight key trends. Omnipod Discover allows users to share data with their Healthcare Professional from using a Clinic code so that their data can be accessed at any time, even between scheduled visits. This means your Healthcare Professional can access your data between visits, providing users are connecting to Wi-fi regularly.[†]

Recommendation: Set up the Controller's Wi-Fi to automatically connect to your home or work network. This will help ensure it connects regularly—ideally as often as possible, but at least once per week.[‡]



During account creation you will be provided with the option to enrol in Omnipod Discover to view and share your Omnipod data.

Click 'Enrol in Discover' if you wish to proceed. Otherwise, you can still continue and complete Omnipod 5 account creation without enrolling to Discover but you must select 'Continue without enrolling' and do not exit the page.



To use Omnipod Discover and share your data you need to click 'open' to complete the following consents:

- Discover Terms & Conditions
- Discover Data Privacy
- Discover Product Research, Development and Improvement Consent



[†]Connections between Omnipod 5 Controller and Insulet Cloud are via Wi-Fi.

To share the data with the medical team and facilitate follow-up, enter the Clinic Code of the Clinic or Healthcare Team.

Share data with your clinic

The ongoing care enter the clinic code to share your diabetes data. The clinic code is a unique identifier to assist you in linking your clinic and uploading sharing.

Note: If you don't know the code you can ask your healthcare provider and add the later on a later online (Omnipod®) account.

Clinic code (see below)

Submit



Share data with your clinic

The ongoing care enter the clinic code to share your diabetes data. The clinic code is a unique identifier to assist you in linking your clinic and uploading sharing.

Note: If you don't know the code you can ask your healthcare provider and add the later on a later online (Omnipod®) account.

Clinic code (see above)

Submit

Clinic group name

Primary clinic address

Agree and share



You are enrolled in Discover!

You have successfully linked your Omnipod® ID to your Discover account and enabled Omnipod® ID System data sharing with this group.

Continue



**OMNIPOD
ACCOUNT
CREATION IS
COMPLETE**

Onboarding Completed!

Omnipod® ID onboarding completed. Your next set up steps will be with your trainer.

Note: Your Omnipod® ID is email address.

You will need your login credentials (Omnipod® ID and Password) to be able to successfully log in to your Omnipod® ID Controller.

Click continue to access your online Omnipod® account.

Continue

Account Linking

Diabetes data partner(s)

Requires your agreement with your Omnipod® ID to allow data to be shared between Omnipod® ID and your Omnipod® ID Controller.

Omnipod® ID System Linking Status: Not linked. To be linked to this controller, please select the Omnipod® ID System you wish to link.

Omnipod® ID System

Link Account



If you decided not to enrol in Omnipod Discover or have not yet received your Clinic code, you can 'Enrol in Discover' at a later date on the Account linking page in your online Omnipod Account.

Questions and Answers

What should I do if I forget my Omnipod ID or password?

The Omnipod ID can be found in the e-mail you received after registration.

If the Omnipod ID or password can't be found, click on "Forgotten ID or password" and follow the steps.



What should I do if I have problems connecting to Omnipod Discover?

If you are experiencing any issues connecting with Omnipod Discover, call 'Customer care'.

What should I do if I haven't received the e-mail?

Check that the e-mail isn't in your spam or junk folders.

Check there is no typo in the e-mail address.

For assistance with account registration please contact Omnipod Customer Care^o:

Kuwait



+965-2461-0486



+965-180-2722



Omnipod-KW@insulet.com



Patient Notes:

Omnipod® ID:

Password:

Clinic Code:

The Omnipod 5 Automated Insulin Delivery System is a single hormone insulin delivery system intended to deliver U-100 insulin subcutaneously for the management of type 1 diabetes in persons aged 2 and older requiring insulin. The Omnipod 5 System is intended for single patient use. The Omnipod 5 System is indicated for use with NovoLog®/NovoRapid®, Humalog®, Trurapi®/Truvelog®/Insulin aspart Sanofi®, Kirsty®, and Admelog/Insulin lispro Sanofi U-100 insulin. Automated mode requires compatible sensor. Compatible sensors are sold and prescribed separately. Sensor availability may vary per country or region.

Refer to the Omnipod® 5 Automated Insulin Delivery System User Guide and www.omnipod.com/safety for complete safety information including indications, contraindications, warnings, cautions, and instructions.

Refer to the Omnipod Discover Data Management System User Guide for complete safety information including

©2026 Insulet Corporation. Insulet, Omnipod, the Omnipod logo, Omnipod Discover and SmartAdjust, are trademarks or registered trademarks of Insulet Corporation. All rights reserved. Dexcom, Dexcom G6, and Dexcom G7 are registered or unregistered trademarks of Dexcom, Inc. in the United States and/or other countries and are used with permission. The sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott and used with permission. The Bluetooth® word marks and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by Insulet Corporation is under license. All other trademarks are the property of their respective owners. The use of third party trademarks does not constitute an endorsement or imply a relationship or other affiliation.

INS-OHS-02-2026-00087 v1.0