

## Multicentre Trial of a Tubeless, On-Body Automated Insulin Delivery System With Customizable Glycemic Targets in Pediatric and Adult Participants With Type 1 Diabetes

- **Clinical objective:** to evaluate the safety and efficacy of the Omnipod<sup>®</sup> 5 Automated Insulin Delivery (AID) System, the first tubeless, on-body AID system with customizable glycemic targets.
- **Primary endpoints** were changed in HbA1c at the end of the AID phase compared with baseline and Time in Range (TIR) 3.9–10.0 mmol/L during the AID phase compared with the Standard Therapy (ST) phase. Primary safety outcomes were incidence of severe hypoglycemia and diabetic ketoacidosis (DKA).
- **Secondary endpoints** included percent time with glucose levels <3.9 mmol/L, >10.0 mmol/L, and <3.0 mmol/L during the AID phase compared with the ST phase.
- **Significant improvements** in HbA1c and glycemic measures, with a low rate of hypoglycemia in a heterogeneous participant group with varied age, baseline glycemia, and prior insulin delivery regimen.

### Study Design


- Prospective, multicentre, single-arm outpatient study:
  - 14-day ST phase
  - 3-month AID phase with Omnipod 5 system
- 240 children, adolescents, and adults enrolled at 17 institutions across the US
- User-selected target glucose ranges from 6.1–8.3 mmol/L
- Unrestricted diet and exercise throughout

### Study Participants

- 112 Children: Age 6 to <14 years
- 128 Adolescents and Adults: Age 14 to 70 years
- All participants:
  - Type 1 diabetes for ≥6 months
  - HbA1c <10.0% (<86 mmol/mol)
  - Any prior insulin therapy (MDI or CSII)

## Omnipod 5 System reduced HbA1c

**HbA1c is reduced by 0.38% (4.2 mmol/mol) in Adolescents and Adults<sup>1</sup>**



Participant endorsement of the system was evident, with 95% enrolling in the extension phase.

● ST phase ● Omnipod 5 System phase

\* $p < 0.0001$ .


<sup>†</sup> $p = 0.046$ .

Baseline and follow-up data were used for the HbA1c primary endpoint. Data shown for Standard Therapy phase and AID phase.

Data shown as mean HbA1c.

MDI, multiple daily injections with insulin.

**HbA1c is reduced by 0.71% (7.8 mmol/mol) in Children<sup>1</sup>**



Connectivity of the on-body devices was excellent, allowing use of automated insulin delivery for median 96.4% of possible time for children.

● ST phase ● Omnipod 5 System phase

\* $p < 0.0001$ .

<sup>†</sup> $p < 0.0005$ .



Baseline and follow-up data were used for the HbA1c primary endpoint. Data shown for Standard Therapy phase and AID phase.

Data shown as mean HbA1c.

MDI, multiple daily injections with insulin.

## Omnipod 5 System increased TIR

- TIR is improved by 2.2 hours/day (9.3%) in Adolescents and Adults<sup>1</sup>
- TIR is improved by 3.7 hours/day (15.6%) in Children<sup>1</sup>



Children spent median 96.4% and adults/adolescents spent median 96.7% of total study time in automated mode.



● ST phase ● Omnipod 5 System phase

\* $p < 0.0001$ .

Data shown as mean % TIR.

## Time spent in hypoglycemia was reduced in adolescents and adults while remaining low in children

- Time in hypoglycemia was reduced by 46% in Adolescents and Adults<sup>1</sup>
- Time in hypoglycemia remained low in Children



● ST phase ● Omnipod 5 System phase

\* $p < 0.0001$ .

Data shown as median % Time Below Ranges <3.0 mmol/L and <3.9 mmol/L.

## Study Highlights<sup>1</sup>:

### Reduced HbA1c

- Omnipod 5 System lowered HbA1c by 0.38% (4.2 mmol/mol) in Adolescents/Adults and by 0.71% (7.8 mmol/mol) in Children
- 66% of Adolescents/Adults and 53% of Children achieved the ADA recommended HbA1c target of <7.0% (<53 mmol/mol)

### Increased TIR

- Improved TIR by 2.2 hours/day (9.3%) in Adolescents/Adults and by 3.7 hours/day (15.6%) in Children
- 69% of Adolescents/Adults and 82% of Children met consensus clinical targets for TIR
- All age groups demonstrated 78% TIR overnight (00:00–06:00 h) during the study

### Low time in hypoglycemia

- Omnipod 5 System reduced time in hypoglycemia by 46% in Adolescents/Adults, including a 60% reduction in nocturnal hypoglycemia (00:00–06:00 h)
- Time in hypoglycemia remained low throughout the study in Children

### System use

- Adolescents/Adults spent a median of 96.7% and Children spent a median of 96.4% of total study time in automated mode
- 3 cases of severe hypoglycemia and 1 case of diabetic ketoacidosis (DKA) were reported in participants during Omnipod 5 System use. These cases were not related to automated insulin delivery malfunction
- The incidence of severe hypoglycemia and DKA were below reported rates in the US T1D Exchange Registry



**Scan code to view  
published study**



This summary has been provided as part of the Omnipod Academy, an educational service provided for Healthcare Professionals by Insulet.

References 1. Adapted from; Brown SA, et al. Multicenter Trial of a Tubeless, On-Body Automated Insulin Delivery System With Customizable Glycemic Targets in Pediatric and Adult Participants With Type 1 Diabetes. *Diabetes Care* 2022; 44:1630-1640.