

# Omnipod Discover™



**Omnipod Discover**

## Provider and Admin Guide

## Contacts and Important Information

### Omnipod Discover Customer Care

<b>Contact Us:</b>	1-800-591-3455 <a href="http://www.omnipod.com/contact-us">www.omnipod.com/contact-us</a>
<b>Address:</b>	Insulet Corporation 100 Nagog Park, Acton MA 01720
<b>Emergency Services:</b>	Dial 911 (USA only; not available in all communities)
<b>Website:</b>	<a href="http://discover.omnipod.com/hcp">discover.omnipod.com/hcp</a>

# Contents

<b>Getting Started with Omnipod Discover</b> .....	5
Indications for Use.....	5
Important Safety Information.....	5
Accessing Omnipod Discover through Omnipod HCP Account.....	6
Register the First Omnipod Discover Clinic Admin.....	6
<b>Standard Staff and Admin Staff Permissions</b> .....	7
<b>Understanding Admin Staff Enrollment Options</b> .....	7
<b>Enroll in Omnipod Discover</b> .....	7
<b>Managing Staff</b> .....	8
Change Permissions.....	8
Invite Staff.....	8
Reinvite Staff.....	9
Approve Staff.....	9
<b>Managing Patients</b> .....	10
How Patient Information is Shared with Your Clinic.....	10
Invite a Patient.....	10
Reinvite a Patient.....	10
Approve Patients.....	11
<b>Uploading Data</b> .....	11
Uploading Data from the Omnipod DASH System to a PC.....	11
Uploading Data from the Omnipod DASH System to a MAC.....	12
<b>About Patient Reports</b> .....	12
Viewing, Printing, and Downloading Reports.....	12
Seeing More than One Date Range.....	13
<b>Omnipod 5 System Report</b> .....	14
Summary.....	14
Daily Report.....	16
Settings.....	18
<b>Omnipod DASH System Report</b> .....	19
Summary.....	19
Daily Report.....	21
Settings.....	23

<b>Glossary</b> .....	<b>24</b>
What Do the Terms on Reports Mean? .....	24
What Do Daily Report Icons Represent? .....	27
Omnipod Discover Label Symbols .....	31
<b>Troubleshooting</b> .....	<b>31</b>
Frequently Asked Questions.....	31
Contact Omnipod Discover Customer Care .....	31
Terms of Use.....	32
HIPAA Privacy Notice .....	32
About Omnipod Discover.....	32

## Getting Started with Omnipod Discover

Omnipod Discover gives you, the healthcare providers (HCPs) and other clinic staff, the ability to view, print and download reports about your Omnipod® 5 or Omnipod DASH® patients' glucose, insulin delivery data and system settings. To get started with Omnipod Discover, a clinic staff member can either respond to an email invitation or log in to their Omnipod HCP Account at [omnipod.com/HCP](https://omnipod.com/HCP) and enroll from there.

Omnipod Discover includes one web application for use by the person with diabetes and another web application for the clinic staff. With Omnipod Discover, the person with diabetes can view weekly data summaries, pattern insights, and/or download the report viewable by their healthcare provider.

### Indications for Use

Omnipod Discover is a retrospective data analytics and reporting system intended for the Omnipod 5 and Omnipod DASH System users or their caregivers and their healthcare providers for the analysis of glucose and insulin delivery data in home and healthcare settings. It is intended as supplemental data for the users to support diabetes management and aid healthcare providers in patient care.

Omnipod Discover is not intended for people with diabetes in acute care settings or for real-time patient monitoring. The Omnipod Discover software platform is not intended to replace the primary real-time display of the sensor or insulin delivery data on the device, nor does it control any functions of the Omnipod system.

Any medical treatment decision or adjustments should not be made based on this software platform; a qualified healthcare professional needs to be consulted to make such decisions.

### Important Safety Information

**Caution:** Omnipod Discover does not display data in real-time and may not include relevant insulin and glucose information needed for real-time treatment decisions. Refer to the Omnipod 5 App or the patient's PDM when making treatment decisions.

**Caution:** Healthcare providers should use information in the software in conjunction with other clinical information available to them.

## Accessing Omnipod Discover through Omnipod HCP Account

Omnipod HCP Account, a web-based portal, provides training materials and more information about Omnipod products.

You can request access to Omnipod Discover by logging into your Omnipod HCP Account at [omnipod.com/hcp](https://omnipod.com/hcp) and clicking “**Enroll in Discover**”. After you have enrolled in Omnipod Discover and the Omnipod Discover admin at your clinic has approved you, then you will see “Go to Discover” when you log in to your Omnipod HCP Account.

### System requirements

You access Omnipod Discover through the Internet using Google Chrome version 95 or later.

### Password requirements

- At least 8 characters.
- An uppercase letter.
- No parts of your username.
- A lowercase letter.
- A number.
- Your password cannot be any of your last 4 passwords.

## Register the First Omnipod Discover Clinic Admin

A clinic needs an Omnipod Discover Admin Staff member before their healthcare providers can use Omnipod Discover.

To register to become your clinic’s first Omnipod Discover Admin Staff member, do one of the following:

- Contact your Omnipod field representative to request to be the first Admin Staff and then follow the instructions from the email you receive.
- You can self-enroll in Omnipod Discover through your Omnipod HCP Account before anyone else from your clinic has enrolled. The first staff member enrolled will automatically be given Admin Staff permissions.

All staff who wish to use Omnipod Discover must have an Omnipod HCP Account. Visit [omnipod.com/hcp](https://omnipod.com/hcp) and follow the on-screen instructions to register for an Omnipod HCP Account.

## Standard Staff and Admin Staff Permissions

Each Omnipod Discover clinic user is given access through permissions. You can have the standard staff permissions or the admin staff permissions. The following table shows what each lets you do.

Permission	Standard Staff	Admin Staff
View/Print/Download Patient Reports.	YES	YES
View Approved Patient Directory.	YES	YES
View Approved Staff Directory.	NO	YES
Approve or Deny Pending Staff.	NO	YES
Approve or Deny Pending Patients.	NO	YES
Invite Patients.	YES	YES
Invite Standard Staff.	NO	YES
Change Staff Permissions.	NO	YES
Remove Staff Access.	NO	YES

## Understanding Admin Staff Enrollment Options

If you have been invited as an Admin Staff from an Insulet Field team member or if you are the first Staff from your clinic to self-enroll, then you will have the Admin Staff permissions.

### Adding more Admin Staff

As an Admin Staff, you can change a Standard Staff member's role to an Admin Staff role.

### Omnipod HCP Account Required

Both Standard Staff and Admin Staff must have an Omnipod HCP Account in order to access Omnipod Discover.

You may be invited via email to create an Omnipod HCP Account or you can go online to create your Omnipod HCP Account.

## Enroll in Omnipod Discover

1. When logged in to your Omnipod HCP Account, click **Enroll Now**.
2. Confirm your clinic email and your clinic name and address, then click **Submit**.
3. Be sure to read the End User License Agreement and check the box if you agree.

**Note:** It may take up to 2 business days to process your enrollment. If you don't receive a registration email within 2 business days, check your spam or junk email folders.

## Managing Staff

### Change Permissions

As an Admin Staff member, you can change permissions of other staff members to Admin Staff or to Standard Staff as follows:

1. Log in to Omnipod Discover with your Admin Staff account.
2. From the Staff directory, click the **Permissions** drop-down for the staff member whose permissions you want to change.
3. Select Admin or Standard.
4. Click **Change** to confirm.

### Enroll or Change Your Omnipod Discover Clinic Admin

Enroll or change your clinic admin as follows:

- Contact Omnipod Product Support or contact your clinic's Omnipod field representative to register as your Omnipod Discover clinic admin or to change your clinic admin. Then see "Register the First Omnipod Discover Clinic Admin" on page 6.

### Invite Staff

You can invite a clinic staff member to join as an Omnipod Discover standard staff member as follows:

1. From Omnipod Discover, select **Staff** under Directory.
2. Click **Invite Staff**.
3. Enter the following for the staff member you want to invite:
  - First Name.
  - Last Name.
  - Professional Email Address.
4. Click **Invite**.

The staff member is added to the directory with the status of "Invited".

## Reinvite Staff

If a staff member's invitation has not yet expired (invitations expire after 14 days), you can reinvite as follows:

1. From the Omnipod Discover Staff Directory, click **Reinvite** in the Status column beside the staff member you want to reinvite.
2. Review to make sure the information is correct.
3. Select **Reinvite**.

**Note:** The "Reinvite" button is only available for invitations that have not expired.

**Tip:** You can also click **Remove Invitation** and then confirm if you no longer want that person to join the clinic.

## Approve Staff

As an Omnipod Discover clinic admin, you can add a staff member who has enrolled in Omnipod Discover through their Omnipod HCP Account to the approved staff directory as follows:

1. Go to [discover.omnipod.com/hcp](https://discover.omnipod.com/hcp). Then, sign in using your Omnipod Discover clinic admin account email address and password.

You see the enrolled Patient list.

2. From the Pending Staff directory, select the staff member(s) you want, then click Approve.

**Note:** If an unknown provider is on the pending list, you may click **Deny**. If you try to invite a staff member who was previously denied, you will be notified of that.

# Managing Patients

## How Patient Information is Shared with Your Clinic

Omnipod 5 and Omnipod DASH users (patients) may give a clinic access to their Omnipod Discover data through onboarding or PodderCentral®. After the patient grants permission to share data with the clinic and a clinic staff member approves, the clinic will be able to view the patient's glucose and insulin data.

## Invite a Patient

Before you begin, enter the patient's last name in the search bar to check if the patient is already in the system.

1. From Omnipod Discover, select **Invite Patient**.

Enter the following for the patient you want to invite:

- First Name.
  - Last Name.
  - Date of birth.
  - Email.
2. Select the patient's Omnipod device.
  3. Select **Invite**.  
The patient is added to the directory with the status of "Invited".

## Reinvite a Patient

If a patient's invitation has not yet expired (invitations expire after 14 days), you can reinvite as follows:

1. From the Omnipod Discover Patient Directory, enter the patient's last name in the search bar.
2. Click **Reinvite** in the Last Updated column beside the patient you want to reinvite.
3. Review to make sure the information is correct.
4. Click **Reinvite**.

**Note:** The "Reinvite" button is only available for invitations that have not expired.

**Tip:** You can also click **Remove Invitation** and then confirm if you no longer want that person to join the clinic.

## Approve Patients

You can approve a pending request to share a patient's information with your clinic as follows:

1. Log in to your Omnipod HCP Account, enter your email address and password, then click **Sign In**.
2. Click **Go to Discover**.
3. From the PENDING REQUESTS section, go to Patients.
4. Select the patient(s) you want, then click **Approve**.

The patient is moved from the Pending list to the Directory.

**Note:** If you don't want to approve, you can select that patient and then click **Deny**.

## Uploading Data

The Omnipod 5 System automatically sends data to Omnipod Discover via Wi-Fi or the embedded SIM card. There is no option to manually upload data from the Omnipod 5 System.

### Omnipod DASH System

Data from the Omnipod DASH System will be uploaded to Omnipod Discover by connecting the PDM to a computer via USB cable and using the upload tool in Omnipod Discover. After upload, you will be able to generate the patient's report. Because the Omnipod DASH System does not include any sensor data, patients will not receive weekly insights. Rather, they will have access to their data and reports at the time of clinic upload only.

**Note:** If the system notices the PDM date and time may be incorrect, it will ask you to enter the date and time shown on the PDM to ensure accurate report generation.

### Uploading Data from the Omnipod DASH System to a PC

1. Connect the PDM to the computer with a cable and export data from the PDM.
2. From Omnipod Discover, search for and select the patient you want to upload data for, then click **Upload Data** for that patient.
3. Confirm you have selected the correct patient, print or read the instructions on how to locate patient data, then click **Upload Patient Data**.
4. In the left panel, find This PC and click **Android Device**.
5. Click **Internal Storage**.
6. Select the OmniPodDash.ibf file and click **Open**.

## Uploading Data from the Omnipod DASH System to a MAC

1. Connect the PDM to the computer with a cable and export data from the PDM. Tap Export on the PDM's screen.
2. Copy and paste the OmniPodDash.ibf patient data file from the Android File Transfer window to your local desktop or other preferred location on your MAC.
3. From Omnipod Discover, search for and select the patient you want to upload data for, then click **Upload Data** for that patient.
4. Confirm you have selected the correct patient, print or read the instructions on how to locate patient data, then click **Upload Patient Data**.
5. Select the OmniPodDash.ibf file from the location you saved it, then open it.

## About Patient Reports

With Omnipod Discover reports, you get a detailed view of how your patient is doing and where they may need additional support.

**Note:** When viewing graphs on your report, look closely at the bolus bars so you can see where multiple boluses are close in time. Seeing bolus instances without a corresponding value attached means that to improve readability, in areas where numbers would be too close so as to overlap, the system removed the numbers where the bolus overlaps.

The report includes the following pages: Summary, Daily Report, and Settings.

## Viewing, Printing, and Downloading Reports

1. Under Directory, in the menu on the left, click **Patients**.
2. Select the patient whose report you want to view.
3. Click the **Select Report Length** dropdown to show the date range options, and select the date range that you want (7,14, 30, 60, or 90 days).

Omnipod Discover will display date range options if valid data is available during those date ranges. For example, if there are 14 days of Omnipod data, the Select Report Length options will offer 7 or 14 day reports, but not 30 days or greater.

**Note:** If there is not enough data available to generate a specific report, you won't be able to select that date range. With less than 7 days worth of data, you will see only a <7 date range option and when viewing reports, certain glucose metrics will not be calculated.

4. Do one of the following:

- Click **View Report** to view the report as a PDF.
- Click **Print Report** if you want to both print the report through your print options window and then view the report as a PDF.
- Click **Download Report** to download the report to your computer.

**Note:** The report file is downloaded to the public download folder for your browser with the report generation date, patient last name, patient first initial, and year of birth as the filename. If you have concerns about possible loss of confidentiality within that folder, move the report out of that folder.

Both **View Report** and **Print Report** options show you the report in PDF format in a new browser tab.

## Seeing More than One Date Range

For Omnipod 5, if you have fewer days of sensor data than your selected report length, you'll see a banner message at the top of the Summary page explaining that glucose-related metrics are based on the connected sensor duration. The affected sections are: Time in Ranges, Glucose Metrics, and Ambulatory Glucose Profile (AGP).

## Signing Out

- Click your name to show the dropdown, then click **Sign Out**.

## Summary - Omnipod 5 System

The Summary section includes both glucose and insulin statistical data, integrated and visualized for your review.

### Time in Ranges

Duration for each range in both percentage and minutes/hours.

### Glucose Metrics & Device Use

Average Glucose

Glucose Management Indicator (GMI)

Glucose Variability (Coefficient of variation)

Standard Deviation (SD)

System Mode Details: % time spent in Automated Mode, Automated Mode: Limited, user-initiated

Activity feature, and Manual Mode

Active Pod Wear Time (%)

Average Pod Wear Time (days/hours)

Connected CGM Time (%)

### Ambulatory Glucose Profile (AGP)

Visualization of hourly glucose variability across the reporting period

### Boluses

See the average number of boluses per day delivered by the user in 6-hour blocks of time:

- Nighttime (12:00 am–5:59am)
- Morning (6:00 am–11:59am)
- Afternoon (12:00 pm–5:59pm)
- Evening (6:00 pm–11:59pm)

### Omnipod 5 Settings

Latest adjustable settings data available from the Omnipod 5 System (Target Glucose, Correct Above, Correction Factor, and Insulin-to- Carb (IC) Ratio) on the same timeline alongside glucose and insulin usage data.

### Insulin Summary

- Average Total Daily Insulin
- Average Daily Basal and Bolus Insulin Usage
- Bolus type, frequency, and amounts
- Diet: Avg carbs/day, Avg # of carb entries/day

## Mary Smith


DOB: Oct 20, 2000 (23 years old)

### 14-Day Summary

May 1, 2024 - May 14, 2024 | Report Generated: Aug 10, 2024

#### Time in Ranges

Goal: <1%  
<54 mg/dL 54-69 mg/dL  
0.4% Very Low 1% Low




#### Glucose Metrics

Average Glucose: 154 mg/dL  
Glucose Mgmt Indicator (GMI): 6.9% Goal: <7%  
Glucose Variability (CV\*): 34.7% Goal: ≤36%  
Standard Deviation (SD): 50 mg/dL  
\*Coefficient of variation

#### Ambulatory Glucose Profile (AGP)

AGP is a summary of glucose values from the report period, with median (50%) and other percentiles shown as if they occurred in a single day.



Bolus (Avg # of boluses/day: 27.2)

Time Period	Avg # of boluses/day
Nighttime	1.7
Morning	7.4
Afternoon	9.9
Evening	8.2

#### Omnipod 5 Settings

Setting	Value
Target Glucose (mg/dL)	120
Correct Above (mg/dL)	130
Insulin to Carb Ratio (unit:grams)	1:10
Correction Factor (unit:mg/dL)	1:15
Target Glucose (mg/dL)	120
Correct Above (mg/dL)	130
Insulin to Carb Ratio (unit:grams)	1:10
Correction Factor (unit:mg/dL)	1:15

#### Insulin Summary

Avg Total Daily Insulin: 59.3 U

Basal/Day: 48% 32.1 U  
Bolus/Day: 52% 27.2 U

Avg # of boluses/day: 27.2  
User override: 4% of total boluses (1 bolus over 14 days)



## Mary Smith

DOB: Oct 20, 2000 (23 years old)

### 14-Day Summary


May 1, 2024 - May 14, 2024 | Report Generated: Aug 10, 2024

Omnipod Discover™

AGP Report - Developed By HealthPartners Institute d/f/a International Diabetes Center

#### Time in Ranges

Goal: <1%  
<54 mg/dL 54-69 mg/dL 70-180 mg/dL 181-250 mg/dL >250 mg/dL  
0.4% Very Low 1% Low 70% Time in Range 23.6% High 5% Very High



#### Glucose Metrics


Average Glucose: 154 mg/dL  
Glucose Mgmt Indicator (GMI): 6.9% Goal: <7%  
Glucose Variability (CV\*): 34.7% Goal: ≤36%  
Standard Deviation (SD): 50 mg/dL  
\*Coefficient of variation

#### Device Use | Omnipod 5 & Dexcom G7

Automated Mode: 96%  
Activity: 4%  
Limited: 3%  
Manual Mode: 4%  
Active Pod Wear Time: 98%  
Avg Pod Wear Time: 2.5 Days (60 hrs)  
Connected CGM Time: 99%

#### Ambulatory Glucose Profile (AGP)

AGP is a summary of glucose values from the report period, with median (50%) and other percentiles shown as if they occurred in a single day.



Bolus (Avg # of boluses/day: 27.2)

Time Period	Avg # of boluses/day
Nighttime	1.7
Morning	7.4
Afternoon	9.9
Evening	8.2

#### Omnipod 5 Settings

Setting	Value
Target Glucose (mg/dL)	120
Correct Above (mg/dL)	130
Insulin to Carb Ratio (unit:grams)	1:10
Correction Factor (unit:mg/dL)	1:15
Target Glucose (mg/dL)	120
Correct Above (mg/dL)	130
Insulin to Carb Ratio (unit:grams)	1:10
Correction Factor (unit:mg/dL)	1:15

#### Insulin Summary

Avg Total Daily Insulin: 59.3 U

Basal/Day: 48% 32.1 U  
Bolus/Day: 52% 27.2 U

Avg # of boluses/day: 27.2  
User override: 4% of total boluses (1 bolus over 14 days)

#### Bolus Type

Meal (may include correction): 62.5% 21.7 U  
Correction Only: 33.3% 3.7 U  
Manual: 4.2% 1.65 U


#### Diet

Avg carbs/day: 183g  
Avg # of carb entries/day: 3.9

## Daily Report - Omnipod 5 System

The Daily Report shows you daily glucose trends, insulin delivery data, carbohydrate intake, and the Omnipod 5 System mode details.

Wednesday, May 1, 2024



**Glucose:** 24-hour tracing shows you glucose values and gaps (if any) in sensor data. A gap or break shown within the otherwise continuous glucose line means that sensor data wasn't available, for example, during a Pod change or Automated Mode: Limited state or due to time zone changes.

**Bolus:** Details and timing of user-initiated boluses for meals and to bring down high glucose.

**Carbs:** Details and timing of carbohydrates entered into the Omnipod 5 App.

## Modes in Omnipod 5 System

<b>Automated Mode</b>	An insulin delivery method in the Omnipod 5 System that uses insulin delivery history, sensor glucose value, and trend to automatically increase, decrease, and pause delivery of insulin based on current and predicted glucose values using a customizable glucose target.	
<b>A</b>	Automated: Limited	Automated insulin delivery in the Omnipod 5 System used when sensor glucose values are not available. Insulin delivery is based on settings and recent insulin delivery.
	<b>L</b>	
	Activity Feature	A user-selected feature in Automated Mode in the Omnipod 5 System for times when less insulin may be needed, such as exercise. When enabled, insulin delivery is reduced, and the Target Glucose is set to 150 mg/dL.
	<b>a</b>	
<b>Manual Mode</b>	Insulin delivery method in the Omnipod 5 System that delivers insulin amounts according to the basal rates in the Basal Program.	
<b>M</b>		

### Additional Basal Indicators

- Max Delivery reached (Automated Mode only)
- Insulin Paused
- Temp Basal (Manual Mode only)

### Pod Changes

- Started new Pod

For additional information about the Omnipod 5 System Modes, please refer to the Omnipod 5 Instructions for Use.


## Mary Smith

DOB: Oct 20, 2000 (23 years old)


Omnipod Discover™

### 14-Day Daily Report


May 1, 2024 - May 14, 2024 | Report Generated: Aug 10, 2024




Wednesday, May 1, 2024



Thursday, May 2, 2024



Friday, May 3, 2024



# Settings page - Omnipod 5 System

The Settings page shows the current Basal and Bolus settings in the Omnipod 5 App.

## Mary Smith

DOB: Oct 20, 2000 (23 years old)

Omnipod Discover™

### Omnipod 5 Settings

SN: 01020000-000000376 | Report Generated: Aug 10, 2024 | This page reflects the most recent pump settings

#### Basal Settings

Changing Basal Programs only impacts Manual Mode basal delivery.

Max Basal Rate 5 U/hr

#### Basal Programs

	Start Time	End Time	Basal Rate
Basal 1 (Active)	12:00 AM	12:00 PM	1.25 U/hr
	12:00 PM	5:00 PM	1.2 U/hr
	5:00 PM	12:00 AM	1.5 U/hr
			Total Basal 31.5 U/day

Temporary Basal Enabled ON

#### Bolus Settings

When using Automated Mode, the main adjustable setting that directly impacts automated insulin delivery is Target Glucose. The SmartBolus Calculator settings below impact suggested bolus doses in both Automated and Manual Mode.

	Start Time	End Time	Target Glucose	Correct Above
Target Glucose & Correct Above	12:00 AM	5:00 AM	120 mg/dL	120 mg/dL
	5:00 AM	8:00 PM	110 mg/dL	110 mg/dL
	8:00 PM	12:00 AM	120 mg/dL	120 mg/dL

	Start Time	End Time	1 U of Insulin Covers
Insulin to Carb (IC) Ratio	12:00 AM	6:00 AM	10 g of carbs
	6:00 AM	12:00 PM	15 g of carbs
	12:00 PM	6:00 PM	12 g of carbs
	6:00 PM	12:00 AM	15 g of carbs

	Start Time	End Time	1 U of Insulin Lowers Glucose By
Correction Factor	12:00 AM	6:00 AM	60 mg/dL
	6:00 AM	11:30 AM	50 mg/dL
	11:30 AM	6:30 PM	40 mg/dL
	6:30 PM	12:00 AM	60 mg/dL

Duration of Insulin Action	4 hours
----------------------------	---------

Max Bolus	25 U
Extended Bolus	ON
Min Glucose for Calculation	70 mg/dL
Reverse Correction	ON

For more information about settings visit [omnipod.com/hcp/clinical-resources/provider-resources](https://omnipod.com/hcp/clinical-resources/provider-resources)

# Summary - Omnipod DASH System

## Glucose Statistics

- Average Glucose
- Number of Glucose Readings
- Lowest - Highest Glucose
- Average Entered Readings per Day
- Glucose Variability (Coefficient of variation)

## Device Use

- Active Pod Wear Time
- Average Pod Wear Time

## Boluses

See the average number of boluses per day delivered by the user in 6-hour blocks of time.

- Nighttime (12:00am - 5:59am)
- Morning (6:00am - 11:59am)
- Afternoon (12:00pm - 5:59pm)
- Evening (6:00pm - 11:59pm)

## Omnipod DASH Settings

Latest adjustable settings data available from the Omnipod DASH System (Target BG, Correct Above, Correction Factor, and Insulin-to-Carb (IC)) on the same timeline alongside Active Basal Program insulin usage.

## Insulin Summary

- Average Total Daily Insulin
- Basal to Bolus (%)
- Bolus: Frequency, amounts, and types

## Diet

- Avg carbs/day
- Avg # of carb entries/day
- Carb presets used

# Sam Smith

DOB: Feb 3, 1995 (28 years old)

## 14-Day Summary

Jul 27, 2024 - Aug 10, 2024 | Report Generated: Aug 10, 2024


### Glucose Statistics

Average Glucose*	154 mg/dL
<small>*Based on entered readings</small>	
Lowest-Highest Glucose	117-243 mg/dL
Glucose Variability (CV*)	54.9% Goal: ≤36%
<small>*Coefficient of variation</small>	


### Device Use | Omnipod DASH

Active Pod Wear Time  
98%

### Bolus



### Omnipod DASH Settings




### Insulin Summary

Avg Total Daily Insulin: 48.1 U



Avg # of boluses/day: 5.4  
User override: 3% of total boluses (4 boluses over 14 days)

### Bolus Type



Omnipod Discover™

# Sam Smith

DOB: Feb 3, 1995 (28 years old)

## 14-Day Summary

Jul 27, 2024 - Aug 10, 2024 | Report Generated: Aug 10, 2024

### Glucose Statistics


Average Glucose*	154 mg/dL	Number of Entered Glucose Readings	39
<small>*Based on entered readings</small>			
Lowest-Highest Glucose	117-243 mg/dL	Average Entered Readings/Day	2.8
Glucose Variability (CV*)	54.9% Goal: ≤36%		
<small>*Coefficient of variation</small>			

### Device Use | Omnipod DASH


Active Pod Wear Time  
98%

Avg Pod Wear Time  
3 Days (72 hours)

### Bolus



### Omnipod DASH Settings




### Insulin Summary

Avg Total Daily Insulin: 48.1 U



Avg # of boluses/day: 5.4  
User override: 3% of total boluses (4 boluses over 14 days)

### Bolus Type



### Diet

Avg carbs/day  
183g


Avg # of carb entries/day  
3.9

Carb presets used  
Yes

## Daily Report - Omnipod DASH System

The Daily Report shows you daily trends in glycemia, insulin delivery, and carbohydrate intake so you can analyze them.

Saturday Jul 27, 2024






**Glucose:** Blood glucose readings from BG meter.


**Bolus:** Details and timing of user-initiated boluses for meals and to bring down high glucose.

**Carbs:** Details and timing of carbohydrates entered into the PDM.

### Basal Activity

-  Scheduled Basal (Omnipod DASH System)
-  Temp Basal (Omnipod DASH System)
-  Suspended (Omnipod DASH System)

### Pod

-  Pod Changed: Started a new Pod

## Sam Smith

DOB: Feb 3, 1995 (28 years old)














Omnipod Discover™

### 14-Day Daily Report


Jul 27, 2024 - Aug 10, 2024 | Report Generated: Aug 10, 2024

Bolus Calculator ON


#### Legend

Glucose Value	Bolus	Basal	Pod
 Above Range	 User Override Up	 Meal Bolus	 Scheduled
 In Target Range	 User Override Down	 Correction Bolus Only	 Temp Basal
 Below Range	 Extended	 Manual Bolus	 Suspend
			 Pod Change


Saturday Jul 27, 2024



Sunday, Jul 28, 2024



Monday, Jul 29, 2024



# Settings page - Omnipod DASH System

The Settings page shows the current Basal and Bolus settings in the Omnipod DASH System.

## Sam Smith

DOB: Feb 3, 1995 (28 years old)

Omnipod Discover™

### Omnipod DASH Settings

SN: 01020000-000000376 | Report Generated: Aug 10, 2024 | This page reflects the most recent pump settings

#### Basal Settings

Max Basal Rate 3 U/hr

##### Basal Programs

	Start Time	End Time	Basal Rate
Basal 1 (Active)	12:00 AM	12:00 PM	1.25 U/hr
	12:00 PM	5:00 PM	1.2 U/hr
	5:00 PM	12:00 AM	1.5 U/hr
			Total Basal 31.5 U/day

Temporary Basal ON

#### Blood Glucose (BG)

	Lower Limit	Upper Limit
BG Goal Limit	70 mg/dL	180 mg/dL

#### Bolus Settings

Bolus Calculator ON

	Start Time	End Time	Target Glucose	Correct Above
Target BG	12:00 AM	5:00 AM	120 mg/dL	120 mg/dL
	5:00 AM	8:00 PM	110 mg/dL	110 mg/dL
	8:00 PM	12:00 AM	120 mg/dL	120 mg/dL

Minimum BG for Bolus Calcs 60 mg/dL

	Start Time	End Time	1 U of Insulin Covers
Insulin to Carb (IC) Ratio	12:00 AM	6:00 AM	10 g of carbs
	6:00 AM	12:00 PM	15 g of carbs
	12:00 PM	6:00 PM	12 g of carbs
	6:00 PM	12:00 AM	15 g of carbs

	Start Time	End Time	1 U of Insulin Lowers Glucose By
Correction Factor	12:00 AM	6:00 AM	60 mg/dL
	6:00 AM	11:30 AM	50 mg/dL
	11:30 AM	6:30 PM	40 mg/dL
	6:30 PM	12:00 AM	60 mg/dL

Reverse Correction ON

Duration of Insulin Action	4 hours
----------------------------	---------

Max Bolus	25 U
Extended Bolus	ON

Carbs	
Carb Presets Used	40g 30g

For more information about settings visit [omnipod.com/hcp/clinical-resources/provider-resources](https://omnipod.com/hcp/clinical-resources/provider-resources)

## Glossary






### What Do the Terms on Reports Mean?






Term	Description
Active Pod Wear Time (%)	% of time the Pod was activated during the report period.
Activity feature	A user-initiated feature in Automated Mode in the Omnipod 5 System for times when less insulin is needed, such as exercise. When enabled, insulin delivery is reduced, and Target Glucose is set to 150 mg/dL.
Afternoon	noon - 5:59pm
Ambulatory Glucose Profile (AGP)	Visualization of hourly glucose variability across the reporting period.
Automated Mode	An insulin delivery method in the Omnipod 5 System that uses insulin delivery history, sensor glucose value, and trend to automatically increase, decrease, and pause delivery of insulin based on current and predicted glucose values using a customizable glucose target.
Automated Mode: Limited	Automated insulin delivery in the Omnipod 5 System used when sensor glucose values are not available. Insulin delivery is based on settings and recent insulin delivery.
Average Entered Readings/Day	Average number of entered glucose readings per day. (Omnipod DASH only)
Average Glucose	Average of all glucose measurements.
Avg Pod Wear Time days/hours	Average number of days/hours each Pod is worn.
Average Total Daily Insulin	Past insulin history that, along with Target Glucose, informs Automated Mode insulin delivery.
Basal to Bolus (%)	How much of total daily insulin goes to basal insulin delivery and how much goes to user-initiated boluses.
Bolus	A dose of insulin delivered for meals with carbohydrates and/or to correct a high glucose.
Bolus Settings	Used to calculate a user-initiated bolus in the bolus calculator.
Carbs (carbohydrates)	Sugars and starches that are consumed and the body breaks down into glucose.
Carb Presets (Custom Foods)	An entry in the Omnipod DASH Food Library's MY FOODS list or the Omnipod 5 Custom Foods list that states the number of carbohydrates contained in a favorite food item, snack, or entire meal. It is assigned a custom name and stored for later use.
Connected CGM Time	% of time sensor (CGM) is connected to the Omnipod 5 Pod.






<b>Term</b>	<b>Description</b>
Correct Above	User setting, where glucose above this lets the bolus calculator calculate a correction bolus.
Correction Bolus	Dose of user-initiated insulin delivered to lower glucose when it is above the target.
Correction Factor	Setting used to indicate how much one unit of insulin will lower blood glucose.
Days in Range	The number of days in a week that time in range was at least 70%.
Duration of Insulin Action	Setting that represents the length of time that insulin stays active in the body. The bolus calculator uses this value to determine Insulin on Board (IOB).
Evening	6:00 pm-11:59 pm
Extended Bolus	Setting (ON or OFF) that provides the ability only during Manual Mode to give a portion of the bolus at the beginning of the meal with the remainder of the bolus dripped over a chosen period of time.
Glucose Management Indicator (GMI)	Provides an estimated A1C level, based on the average glucose level from sensor readings for 14 or more days of data.
Glucose Variability (Coefficient of variation)	The measurement of fluctuations of glucose over time.
Insulin Paused	Insulin delivery stopped.
Insulin to Carb (IC) Ratio	How many carbs covered by one unit of insulin.
Manual Bolus	A bolus amount chosen by the user (not calculated by the bolus calculator).
Manual Mode	Insulin delivery method that delivers insulin amounts according to the basal rates in the Basal Program.
Max Basal Rate	Defines an upper limit for any basal rate used in Basal Programs and temp basals only during Manual Mode.
Max Bolus	Defines the upper limit for a user-initiated bolus. The bolus calculator prevents the user from entering a bolus over this amount. It also prevents the user from entering a carb amount that will cause the calculated bolus to exceed this amount. The largest allowed value for Maximum Bolus is 30 units.
Meal Bolus	Dose of insulin taken to cover carbohydrates in a meal or snack.
Min Glucose for Calculation	When glucose is below this value, the SmartBolus Calculator is disabled and does not calculate a bolus.
Mode	The Omnipod 5 System Mode (Automated, Manual, Automated: Limited, and Activity feature).
Morning	6:00 am - 11:59 am
Nighttime	Midnight - 5:59 am







<b>Term</b>	<b>Description</b>
Number of Entered Glucose Readings	Total number of entered glucose readings for the report period.
Reverse Correction	Setting (ON or OFF) that determines how the bolus calculator handles meal boluses when glucose is below Target Glucose.
Sensor-based glucose monitoring system	System to track glucose throughout the day and night, supplied by a third-party medical device manufacturer. Often called CGM or Sensor.
Sensor glucose value	Glucose measured by a Sensor. Sensor glucose values include trend, which indicates whether glucose is going up, down, or remaining steady.
Settings	Most recent uploaded adjustable settings programmed into the Omnipod 5 Controller on the same timeline alongside glucose and insulin usage data.
SmartAdjust™ technology	Pod software used to calculate automated insulin delivery as often as every 5 minutes to bring the user's glucose to the customized glucose target, or Target Glucose.
System Mode Details	% time in Automated Mode, Automated Mode: Limited, Automated Mode with Activity feature, and Manual Mode.
Target Glucose	The user-customizable glucose target used by both SmartAdjust technology and the Omnipod 5 System SmartBolus Calculator to calculate how much insulin is needed based on both manually entered blood glucose readings and the glucose values from the sensor. Target Glucose can be set from 110-150 mg/dL in 10 mg/dL increments and in up to 8 Target Glucose segments per day.
Temp Basal	Temporary basal rate used to handle a change in routine.
Temporary Basal Enabled	Lets the user override the currently running Basal Program by increasing or decreasing it for a predetermined period of time (If using Omnipod 5, available only in Manual Mode).
Time in Mode	Time in one of these modes: Automated, Manual, Automated: Limited, and Automated with user-initiated Activity feature.
Time in Range	The percent of time with glucose in a target range of 70–180 mg/dL. An increase in time in range means less time spent with high and low glucose. A 4% increase or decrease in time in range is about 1 hour.
User Override	When a user changes the bolus amount suggested by the bolus calculator to deliver more or less bolus insulin.

## What Do Daily Report Icons Represent?






Icons on daily report	What it means
	Sensor Glucose Value and Trend
	Manually Entered Glucose
	Target Glucose
	User Override Up
	User Override Down

Icons on daily report	What it means
	<p>Extended Bolus</p>
	<p>Meal Bolus</p>
	<p>Correction Bolus</p>
	<p>Manual Bolus</p>
	<p><b>Automated Mode:</b> An insulin delivery method in the Omnipod 5 System that uses insulin delivery history, sensor glucose value and trend to automatically increase, decrease, and pause delivery of insulin based on current and predicted glucose values using a customizable glucose target.</p>

Icons on daily report	What it means
	<p><b>Activity Feature:</b> A user-selected feature in Automated Mode in the Omnipod 5 System for times when less insulin may be needed, such as exercise. When enabled, insulin delivery is reduced, and the Target Glucose is set to 150 mg/dL.</p>
	<p><b>Automated Mode: Limited:</b> Automated insulin delivery in the Omnipod 5 System used when sensor glucose values are not available. Insulin delivery is based on user settings and recent insulin delivery.</p>
	<p><b>Manual Mode:</b> Insulin delivery method in the Omnipod 5 System that delivers insulin amounts according to the basal rates in the user's Basal Program.</p>
	<p>Maximum Insulin Delivery (Automated Mode only)</p>
	<p>Insulin Paused (Omnipod 5 System)</p>

Icons on daily report	What it means
	Temp Basal (Omnipod 5 System)
	Insulin Suspended (Omnipod DASH System)
	Scheduled Basal (Omnipod DASH System)
	Temp Basal (Omnipod DASH System)
	Bolus when the bolus calculator is off (Omnipod DASH System).
	<b>Pod Changed:</b> Started a new Pod

## Omnipod Discover Label Symbols

Symbol	Meaning	Symbol	Meaning
	Date of manufacture		Manufacturer
	Catalogue number		Unique device identifier
	Consult instructions for use or consult electronic instructions for use	<b>Rx ONLY</b>	Caution: Federal law restricts this device to sale by or on the order of a physician.

## Troubleshooting

### Frequently Asked Questions

Issue or Error message	What you can do
Can't sign in	Tap <b>Forgot username/password</b> on the "Sign in" screen to reset your password. If you continue to have trouble, contact Customer Care.  Wait at least 1 day after enrollment to sign in.
Patient not on list	Follow instructions in this guide to invite the patient. After the patient enrolls in Omnipod Discover, if you still cannot view their data, contact Customer Care.
Missing data	Check the Wi-Fi or Cellular connectivity settings on the patient's Controller. If you are still experiencing this issue, contact Customer Care.

### Contact Omnipod Discover Customer Care

If you have a problem with using Omnipod Discover, contact Omnipod Discover Customer Care at 1-800-591-3455.

## Terms of Use

The use of the Omnipod Discover platform, including any reports generated by the system as well as any information contained in this guide, is subject to the Omnipod Discover Terms of Use, including the End User License Agreement available at:

<https://discover.omnipod.com/docs/legal.html>

and Cookie Policy available at:

<https://discover.omnipod.com/docs/cookies.html>

## HIPAA Privacy Notice

In addition to the Omnipod Discover Terms of Use, the use of the services is also subject to Insulet's HIPAA Privacy Notice, available at:

<https://www.omnipod.com/hipaa-notice>

## About Omnipod Discover

### Version

1.5



Insulet Corporation

100 Nagog Park, Acton, MA 01720 USA

omnipod.com

# Omnipod Discover™



**Insulet Corporation**

100 Nagog Park  
Acton, MA 01720 USA  
1-800-591-3455  
**omnipod.com**

© 2024 Insulet Corporation. Insulet, Omnipod, the Omnipod logo, Omnipod Discover, DASH, and SmartAdjust are trademarks or registered trademarks of Insulet Corporation. All rights reserved. All other trademarks are the property of their respective owners. The use of third party trademarks does not constitute an endorsement or imply a relationship or other affiliation.

Patent information at [www.insulet.com/patents](http://www.insulet.com/patents)

PT-000944-AW Rev 09 2024-11-21