Omnipod® 5 is now compatible with the Dexcom® G7 receiver and the Dexcom G7 app.

Dexcom has released a new version of the Dexcom G7 sensor, which allows Omnipod 5 users to manage their sensor using either the Dexcom G7 app or the Dexcom G7 receiver (no longer requiring exclusive use of the Dexcom G7 app).

As part of this update, the Dexcom G7 sensor will send glucose values to both the Pod and the Dexcom display device (either the Dexcom G7 app or the Dexcom G7 receiver).

Users should follow the updated instructions below to ensure proper Pod–sensor connectivity:

Updated Instructions for Dexcom G7 users

- The connection principle remains unchanged: the Dexcom G7 Sensor must be connected to both the Omnipod 5 System and a Dexcom display device, which can now be either the Dexcom G7 app or the Dexcom G7 receiver.