

Omnipod® Comparison Chart

	Features	omnipod	omnipod DASH
The Pod	Pod Size (cm)	3.9 W x 5.2 L x 1.45 H	3.9 W x 5.2 L x 1.45 H
	Pod Weight (g without insulin)	25	26
	Reservoir	200 units	200 units
	Infusion Set	Automated cannula insertion	Automated cannula insertion
	Waterproof	Yes*	Yes*
The Personal Diabetes Manager (PDM)	PDM Size (cm)	6.21 W x 11.25 L x 2.5 H	6.3 W x 13.0 L x 1.0 H
	PDM Weight (g)	125	106
	Touch Screen	No	Yes
	Battery Type	Disposable AAA alkaline (2)	Rechargeable lithium ion battery
	Battery Life	Approximately 3 weeks	Approximately 2 days of typical use following a full charge
	CGM Integration	No	No
	BGM Compatibility	Built into PDM	Manual Entry. Any BGM can be used
	Pod Site Tracker	No	Yes
	Data Management Software	Glooko®	Glooko®
	Wireless Source	Bidirectional radio frequency	Bluetooth® Wireless Technology
	Range	<ul style="list-style-type: none"> Up to 1.5m to deliver a bolus Pod and PDM must be side by side to activate a new Pod Not required to be near the Pod when the Pod is active and delivering basal insulin 	<ul style="list-style-type: none"> Up to 1.5m to deliver a bolus Pod and PDM must be side by side to activate a new Pod Not required to be near the Pod when the Pod is active and delivering basal insulin
	Insulin Delivery	Basal Increments	0.05 u/hr
Zero Basal Rate		No	Yes
Insulin-to-Carb Ratio		1 g/unit of insulin	0.1 g/unit of insulin
Basal Programs		7	12
Temp Basal Presets		7	12
Bolus Presets		7	7

* THE POD HAS A WATERPROOF IP28 RATING FOR UP TO 25 FEET FOR 60 MINUTES. THE PERSONAL DIABETES MANAGER IS NOT WATERPROOF.



Insulet Canada Corporation
1540 Cornwall Rd, Suite 201
Oakville, ON, L6J 7W5
1-855-POD-INFO (763-4636)
www.omnipod.com

© 2022 INSULET CORPORATION. OMNIPOD, THE OMNIPOD LOGO, PODDER, DASH, THE DASH LOGO, AND SIMPLIFY LIFE ARE TRADEMARKS OR REGISTERED TRADEMARKS OF INSULET CORPORATION. ALL RIGHTS RESERVED. GLOOKO IS A TRADEMARK OF GLOOKO, INC. AND USED WITH PERMISSION. THE BLUETOOTH® WORDMARK AND LOGOS ARE REGISTERED TRADEMARKS OWNED BY THE BLUETOOTH SIG, INC. AND ANY USE OF SUCH WORDMARKS BY INSULET CORPORATION IS UNDER LICENCE. ALL OTHER TRADEMARKS ARE THE PROPERTY OF THEIR RESPECTIVE OWNERS. THE USE OF THIRD PARTY TRADEMARKS DOES NOT CONSTITUTE AN ENDORSEMENT OR IMPLY A RELATIONSHIP OR OTHER AFFILIATION.