

# THERE ARE LOTS OF PUMPS. THERE'S ONLY ONE OMNIPOD® 5



Screen image is an example, for illustrative purposes only.  
Pod shown without necessary adhesive.



## OMNIPOD® 5 SIMPLIFY LIFE®

The Pod has a waterproof IP28 rating for up to 7.6 metres (25 feet) for 60 minutes.  
The Controller is not waterproof.



### References:

1. Brown S. et al. Diabetes Care. 2021;44:1630-1640. Prospective pivotal trial in 240 participants with T1D aged 6 - 70 yrs. Study included a 14-day standard therapy (ST) phase followed by a 3-month Omnipod 5 hybrid closed-loop (HCL) phase. Mean time in range (3.9-10.0 mmol/L or 70-180mg/dL) in adults/adolescents as measured by CGM: ST = 64.7%, 3-mo Omnipod 5 = 73.9%, P<0.0001.
2. Sherr J. et al. Diabetes Care. 2022; 45:1907-1910. Single-arm multicenter clinical trial in 80 pre-school children (aged 2-5.9 yrs) with T1D. Study included a 14-day standard therapy (ST) phase followed by a 3-month AID phase with Omnipod 5 system. Mean time in range (3.9-10.0 mmol/L or 70-180mg/dL) as measured by CGM in children ST vs. 3-mo Omnipod 5: 57.2% vs 68.1%, P<0.0001.

The Pod has a waterproof IP28 rating for up to 7.6 metres (25 feet) for 60 minutes. The Controller is not waterproof. The Omnipod 5 Automated Insulin Delivery System is a single hormone insulin delivery system intended to deliver U-100 insulin subcutaneously for the management of type 1 diabetes in persons aged 2 and older requiring insulin. The Omnipod 5 System is intended for single patient use. The Omnipod 5 System is indicated for use with NovoLog®/NovoRapid®, Humalog®, Trurapi®/Truvelog®/Insulin aspart Sanofi®, Kirsty®, and Admelog/Insulin lispro Sanofi U-100 insulin.

Refer to the Omnipod® 5 Automated Insulin Delivery System User Guide and [www.omnipod.com/safety](http://www.omnipod.com/safety) for complete safety information including indications, contraindications, warnings, cautions, and instructions. Refer to the Omnipod Discover Data Management System User Guide for complete safety information including indications, contraindications, warnings, cautions, and instructions.

©2025 Insulet Corporation. Omnipod, the Omnipod logo, Simplify Life and SmartAdjust are trademarks or registered trademarks of Insulet Corporation. All rights reserved. The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by Insulet Corporation is under license. All other trademarks are the property of their respective owners. The use of third party trademarks does not constitute an endorsement or imply a relationship or other affiliation. INS-OHS-11-2025-00105 V1.0  
MDMA-2-2024-1385

# MAKE DIABETES MANAGEMENT A SMALLER PART OF LIFE



For people with  
type 1 diabetes,  
ages 2 years  
and older



# TUBELESS, AUTOMATED INSULIN DELIVERY

**Omnipod® 5 is a two part system that integrates with compatible sensors**



## Omnipod 5 Controller

Take control with the Omnipod 5 handheld device.

## Pod

Tubeless, wearable and waterproof\*\* with built-in SmartAdjust™ technology — a predictive insulin delivery algorithm — the Pod sits right on your body and helps to keep glucose in range.<sup>1,2</sup>

## Your sensor

In constant communication, your sensor\* sends glucose values and trends to the Pod to help the System make automated insulin delivery decisions. Please visit [omnipod.com](http://omnipod.com) for a list of compatible sensors.

Screen image is an example, for illustrative purposes only. Pod shown without necessary adhesive.  
\*Compatible sensors are sold and prescribed separately. Sensor availability may vary per country or region. \*\*The Pod has a waterproof IP28 rating for up to 7.6 metres (25 feet) for 60 minutes. The Controller is not waterproof.

# ADJUSTING SO YOU DON'T HAVE TO

**SmartAdjust™ technology** is the brains behind the Omnipod 5 System, helping to protect against highs and lows.<sup>1,2</sup> The System predicts where glucose values will be in 60 minutes and automatically increases, decreases, or pauses insulin every 5 minutes based on your sensor value and trend. Bolus doses can be given for meals using the Omnipod 5 Controller.



# MEALTIME MATHS, DONE FOR YOU

The SmartBolus Calculator makes it easy to determine the right amount of insulin for meals and snacks. It uses your sensor value and trend along with carbohydrate amount to calculate a precise mealtime dose – no need for guessing. Plus, with just a few taps on a touchscreen, it's simple to teach the other caregivers in their life to use – like grandparents, the babysitter, or their school nurse.



## GET GOING WITH ACTIVITY FEATURE

When turned on, Activity Feature reduces insulin delivery for times when glucose typically goes low, like exercising. Whether hitting your step goals, the gym or swimming\*\*, Omnipod 5 makes the balancing act easier.



## CONTROLLED BY A HANDHELD DEVICE

Take full control of the System at any time with just a few taps. Bolusing for a meal, enabling Activity Feature and even activating your Pod looks just like sending a text with the Omnipod 5 Controller.