

Welcome to Omnipod® 5!


New to Omnipod® 5? Enrol in Omnipod Discover™




Omnipod Discover™ is the new, official data tool from Omnipod that delivers weekly, personalised insights into your glucose and insulin patterns to simplify your daily management. By sharing these intuitive reports with your healthcare professional, you can make every check-up easier.

To begin:


Scan this code on your smartphone or visit omnipod.com/setup




Select 'Let's Start'.




Sign in or create your Omnipod ID.




Review and agree to the Terms & Conditions and Omnipod® 5 consents.




Complete the steps to Enrol in Omnipod Discover™ to share diabetes data with your healthcare professional.




Enter the clinic code provided by your healthcare professional.



After enrolling in Omnipod Discover, Tap 'Continue' to access your online Omnipod account.



Review and agree to the Omnipod® Account Access Consent.



Clinic Code:



Already using Omnipod® 5?

Enrol in Omnipod Discover™




To begin:

Scan this code on your smartphone
or visit: account-intl.omnipod.com

Share your diabetes data with your healthcare professionals and access data insights to help you manage your diabetes in just 3 easy steps:

- 1 Sign into your Omnipod account, navigate to **My Account** and select **Account Linking** in the main menu
- 2 Click **‘Enrol in Discover’**
- 3 When asked, **enter the clinic code** to share data with your healthcare professional.

Clinic Code:



Always follow the directions for use.

The Omnipod® 5 Automated Insulin Delivery System is a single hormone insulin delivery system intended to deliver U-100 insulin subcutaneously for the management of type 1 diabetes in persons aged 2 and older requiring insulin.

Refer to the Omnipod® 5 Automated Insulin Delivery System User Guide and www.omnipod.com/en-au/important/safety/information for complete safety information including indications, contraindications, warnings, cautions, and instructions.

Omnipod Discover™ is a retrospective data analytics and reporting system intended for Omnipod® 5 System users or their caregivers and their healthcare providers for the analysis of glucose and insulin delivery data in home and healthcare settings. It is intended as supplemental data for the users to support diabetes management and aid healthcare providers in patient care. Omnipod Discover™ is not intended for people with diabetes in acute care settings or for real-time patient monitoring. The Omnipod Discover software platform is not intended to replace the primary real-time display of the sensor or insulin delivery data on the device, nor does it control any functions of the Omnipod® 5 System. Any medical treatment decision or adjustments should not be made based on this software platform; a qualified healthcare professional needs to be consulted to make such decisions.

©2026 Insulet Corporation. Insulet, Omnipod, the Omnipod logo, are trademarks or registered trademarks of Insulet Corporation. All rights reserved. Insulet Australia Pty Ltd, Level 16, Tower 2, Darling Park, 201 Sussex Street, Sydney, NSW, 2000. 1800 954 075 INS-OHS-03-2026-00170 V 1.0